

METREL Medical ES Manager Instruction manual

Version 1.1.2 © 2023 METREL



Table of Contents

1	Welcome to Metrel Me	dical ES Manager Help	5
2	2 Introduction to Metrel	Medical ES Manager	6
	2.1 What is Metrel Me	edical ES Manager?	6
3	B End-User Licence Agre	eement	7
	-		
	3.1.1 Licence grant		7
	LICENSE AGREEMEN	IT	7
4	Software Installation		10
	4.1 System Requirem	ients	10
	5 1	tware	
	4.3 Metrel Medical ES	S Manager first use	12
	4.4 Software updates	5	12
5	5 Getting started		14
	5.1 Quick guides and	tutorials	14
	5	uides	
6	Liser Interface Compon	nents	15
Ű	•		
		J functions	
		Menu	
	•	tions	
	6.2.2.1 Documen	nt	18
		ication	
		unctions	
		nu group	
		e menu group ments menu group	
		nenu group	
	5	menu	
	0		
	6.2.4.1 Contacts	data lists	23
	6.2.4.2 Structure	e Names data list	24
		_ists	
		S Manager User workspace	
	1	e management	
		DNS	
		rea options /iew filters	
		dow options	

7	Creating Data Structures	
	7.1 Considerations on a new test project's data structure	
	7.2 Creating a new data structure file for a test project	
	7.2.1 Creating a new data structure	
	7.2.2 Appending measurements to the data structure elements	
	7.2.2.1 Appending Single tests	
	7.2.2.2 Appending Auto Sequence [®] test	
	7.3 Creating Data structure using Copy and Paste commands	
	7.3.1 Example of Copy and Paste within the same data structure file	
	7.3.2 Copy and Paste between different data structure files	
	7.3.3 Example of merge functionality	
	7.3.4 Example of copy and paste as same functionality	41
8	Communication with the instruments	
	8.1 Establishing communication with instruments	
	8.1.1 Instrument info	
	8.1.2 Upgrading the Firmware	
	8.2 Upload or download Project files	
	8.2.1 Single file transfer options	
	8.2.2 Multiple file transfer options	
9	Measurements	
	9.1 Single tests	
	9.1.1 Single tests in the user workspace	
	9.1.2 Viewing or configuring the single tests	
	9.1.3 Running single tests	
	9.1.4 Single test results	
	9.2 Auto Sequence [®] tests	
	9.2.1 Auto Sequence [®] tests in the user workspace	
	9.2.2 Viewing and configuring an Auto Sequence [®]	
	9.2.3 Running Auto Sequences [®]	
	9.2.4 Auto Sequence results	
	9.3 Standardized inspections	
	9.3.1 Standardized inspection test in the user workspace	
	9.3.2 Executing the standardized inspection tests	54
10) Measurement Organizer	
	10.1 Measurements Organizer main menu	
	10.2 Auto Sequence [®] Editor Workspace	57
	10.2.1 Auto Sequence [®] group tree structure	57
	10.2.2 Auto Sequence [®] content view	
	10.2.2.1 Auto Sequence [®] Header step options	
	10.2.2.2 Auto Sequence [®] Measurement step	61
	10.2.2.3 Auto Sequence® Result screen step	
	10.2.3 Flow commands	
	10.2.3.1 Number of measurement steps	
	10.2.3.2 Pause	
	10.2.3.3 Buzzer mode	
	10.2.3.4 No notifications mode	
	10.2.3.5 Inspection Expert mode	64
	10.2.3.6 Operation after end of test	

10.3	Custom Inspection Editor Workspace	64
11 Print	ting results	67
11.1	Printing results options	67
11.2	Printing results format	67
12 Crea	ting Reports	70
13 Trou	ıbleshooting	71
13.1	Introduction to troubleshooting	71
13.2	Instrument connection troubleshooting	
13.3	Reporting problems to Metrel	71
13.3.1	1 Bug report checklist	72
14 Onlir	ne support information	73
14.1	How to get support?	73
14.1.1		
14.1.2	2 Use the online contact form	73
Appendix	A Using Auto Sequence®s to test according to IEC/EN 60601	74
Getting	started	74
Measure	ements and parameters	76
Example	e: patient monitor Philips MP70	77

1 Welcome to Metrel Medical ES Manager Help

Use the Bookmarks on the left side of the window to select or search for a topic. You can also select one of the following **Quick links** to start learning about ES Manager.

- ► <u>Introduction</u>
- Online support information
- End-User Licence Agreement
- Getting started



Note

Please note that some parts of program can be changed without notice in new software versions, and therefore might differ from the information in this manual. If you encounter such differences, we will gladly accept your e-mail with your suggestions and include or change the necessary information in our next hotfix or update

2 Introduction to Metrel Medical ES Manager

2.1 What is Metrel Medical ES Manager?

Metrel Medical Electrical Safety Manager is a PC software application for management of safety of medical equipment.

It can be used in combination with new Metrel's test equipment for medical environment. Medical ES Manager enables:

- Pre-treatment of test element structures on the PC.
- Uploading and downloading test results to/from the instrument, downloading.
- Viewing, editing, archiving of test data.
- Generating and printing test reports.
- Creating custom Auto Sequence® tests and custom inspection test.
- Carrying out measurements directly from the PC
- Connection to Metrel Cloud / Webreports

Functionality of the Metrel Medical ES Manager PC software program depends on the licence key, which is always related to the serial number of the purchased instrument.

3 End-User Licence Agreement

3.1 Important Notice

Read this License Agreement carefully before using this Software. BY USING THIS SOFTWARE IN ANY WAY, YOU ACKNOWLEDGE THAT YOU HAVE READ, UNDERSTOOD, AND AGREED TO THE TERMS OF THIS AGREEMENT. IF YOU DO NOT AGREE TO THESE TERMS, DO NOT USE THIS SOFTWARE IN ANY WAY, AND PROMPTLY RETURN IT OR DELETE ANY COPIES OF THIS SOFTWARE IN YOUR POSSESSION.

3.1.1 Licence grant

By installing the software Metrel Medical ES Manager (hereinafter "the Software" or "Software"), you are accepting the following License Agreement.

LICENSE AGREEMENT

This is a legal agreement between You (either an individual or an entity) and Metrel d.o.o. By installing the software or using it in any way you are agreeing to be bound by the terms of this agreement (hereinafter "Agreement"). If you do not agree to the terms of this agreement do not install or use the Software in any way.

- GRANT OF LICENSE. Metrel d.o.o. grants to You (including your employees and subcontractors) a revocable, non-exclusive and limited License to make and use copies of the Software strictly in accordance with the terms and conditions of this Agreement. The rights granted under this Agreement are only effective upon payment of the License fees. The Software is licensed, not sold. The Software, its documentation and other intellectual property related to it are owned by Metrel d.o.o. and protected by the international copyright laws. Metrel grants You the License to:
 - a. Install and run the software on a single or multiple computers used to connect to Metrel instruments. The License is tied to each instrument's serial number. Each instrument requires a separately purchased License.
 - b. Make one copy of the Software solely for backup or archival purposes.
 - c. Transfer the Software to a single hard disk provided you keep the original solely for backup or archival purposes.
 - d. You may store or install a copy of the Software on a storage device, such as a network server, used only to install or run the Software on other computers over an internal network; however, you must purchase and dedicate a License for each separate computer on which the Software is installed or run from the storage device. A License for the Software may not be shared or used concurrently on different computers.
- 2. RESTRICTIONS OF USE. You shall use the Software in accordance with the terms of this Agreement and you shall not:
 - a. Decompile, reverse engineer, disassemble, attempt to derive the source code for or decrypt the Software.
 - b. Make any modification, adaptation, improvement, enhancement, derivative or translation of the Software.
 - c. Violate any applicable law, rule or regulation in connection to your access to the Software.
 - d. Remove, obscure or alter any proprietary notice (including copyright) of Metrel d.o.o. or its partners, except in reports where changing the default logo is specifically enabled by the **Software's functions**.
 - e. Use the Software for any application it is not designed and intended for.

- f. Use the Software to create a product, service or software that is (directly or indirectly) in competition or a substitute for any offered by Metrel d.o.o.
- g. Use the Software to send automated queries to any website or to send any unsolicited commercial messages.
- h. Use any proprietary information, user interfaces or other intellectual property of Metrel d.o.o. in the design, development, manufacture, licencing or distribution of any new software, accessory or device.
- 3. TERMINATION. This Agreement remains effective from Your implicit or explicit agreement to it until termination. It will immediately terminate upon notice if You exceed the scope of License granted or otherwise fail to comply to this Agreement. You may terminate this Agreement by uninstalling the Software, returning all purchased License keys and ceasing the use of the Software and its documentation.
- 4. NO WARRANTY. EXCEPT AS EXPRESSLY SET FORTH IN THIS AGREEMENT, THE PROGRAM IS PROVIDED ON AN "AS IS" BASIS, WITHOUT WARRANTIES OR CONDITIONS OF ANY KIND, EITHER EXPRESS OR IMPLIED INCLUDING, WITHOUT LIMITATION, ANY WARRANTIES OR CONDITIONS OF TITLE, NON-INFRINGEMENT, MERCHANTABILITY OR FITNESS FOR A PARTICULAR PURPOSE. Each recipient is solely responsible for determining the appropriateness of using and distributing the software and assumes all risks associated with its exercise of rights under this Agreement, including but not limited to the risks and costs of program errors, compliance with applicable laws, damage to or loss of data, programs or equipment, and unavailability or interruption of operations.
- 5. DISCLAIMER OF LIABILITY. EXCEPT AS EXPRESSLY SET FORTH IN THIS AGREEMENT, NEITHER RECIPIENT NOR ANY CONTRIBUTORS SHALL HAVE ANY LIABILITY FOR ANY DIRECT, INDIRECT, INCIDENTAL, SPECIAL, EXEMPLARY, OR CONSEQUENTIAL DAMAGES (INCLUDING WITHOUT LIMITATION LOST PROFITS), HOWEVER CAUSED AND ON ANY THEORY OF LIABILITY, WHETHER IN CONTRACT, STRICT LIABILITY, OR TORT (INCLUDING NEGLIGENCE OR OTHERWISE) ARISING IN ANY WAY OUT OF THE USE OR DISTRIBUTION OF THE PROGRAM OR THE EXERCISE OF ANY RIGHTS GRANTED HEREUNDER, EVEN IF ADVISED OF THE POSSIBILITY OF SUCH DAMAGES.
- 6. COPYRIGHT. Any and all trademarks or service marks that Metrel d.o.o. uses in connection with the Software or with services rendered by Metrel d.o.o., are marks owned by Metrel d.o.o.. This Agreement does not grant you any right, license, or interest in such marks, and you shall not assert any right, license, or interest in such marks or designs that are confusingly similar to such marks.
- 7. PRIVACY POLICY.

Metrel d.o.o . will only act upon the Data Controller's instructions with regard to processing data in accordance with GDPR. Metrel d.o.o . will comply with all applicable laws pertaining to personal data and metadata including GDPR. Where requested Metrel d.o.o. are prepared to enter into separate agreements regarding data processing

The Data Controller is defined in the Metrel's privacy policy available on https://www.metrel.si/en/content/privacy/android-privacy-policy.html.

Metrel d.o.o. may (and you hereby expressly authorize Metrel d.o.o. to) disclose information about you to law enforcement officers, as Metrel d.o.o. in its sole discretion, deem necessary or appropriate to investigate or resolve possible crimes or to respond to judicial, regulatory, agency or similar inquiries. Metrel d.o.o. may also contact you using the information you provided to, for example, provide technical

support, respond to user inquiries and transmit questionnaires".

You acknowledge receipt of and agree to Metrel's privacy policy which is made available in full at <u>https://www.metrel.si/en/content/privacy/android-privacy-policy.html</u>. You hereby expressly consent to Metrel's collection, storage and processing of your personal data according to Metrel's current privacy

policy as of the Effective Date, which is incorporated into this End User Agreement by reference. Metrel publishes a privacy policy on its <u>https://www.metrel.si/en/content/privacy/android-privacy-policy.html</u> and may amend such policy from time to time in its sole discretion. You acknowledge that that you had the opportunity to review a more detailed explanation of how your information will be collected, stored and processed by Metrel in accordance with the published privacy policy of Metrel and that you consented to Metrel's privacy policy prior to agreeing to this Agreement. You will ensure that each member of your organization (including employees and contractors) about whom personal data may be provided to Metrel has given his or her express consent to Met**rel's processing of such personal** data.

4 Software Installation

4.1 System Requirements

Before installation, you should check that your system meets the following requirements:

Supported Operating Systems:

- Windows 10, 32-bit or 64-bit
- Windows 11

Installed System Memory (RAM):

• 4 GB (8 GB recommended)

Hard Disk Space:

- At least 650 MB of free space for the installation files and documentation
- Additional 280 MB (x86) or 610 MB (x64) of free space if Microsoft.NET Framework (4.0 or higher) is not previously installed
- Additional disk space (20 GB recommended) is needed for downloading and saving Medical ES Manager data files (A single Metrel Medical ES Manager data file can be up to 8 GB). Additional disk space (20 GB recommended) is needed for downloading and saving Metrel Medical ES Manager data.

4.2 Installing the software

To install the software, download the application from Metrel Download centre and start the installation by running setup.exe.

The installation wizard will guide you through the selection of language used during installation procedure, and installation folder selection (where Metrel Medical ES Manager program will be installed). The default location offered is 'C:\Program Files (x86)\Metrel\Medical MESM ' Browse for your preferred location if you want to change it.

By default the option 'Associate Metrel Medical ES Manager with .padfx file extension' is checked. Metrel Medical ES Manager will be started automatically when data structure file (*filename.padfx*) is opened from PC file manager.

Check 'Create a desktop shortcut' to add an icon for starting Metrel Medical ES Manager from the Windows desktop.

A survey of selected Metrel Medical ES Manager options is presented before the starting the installation. See *Figure 4.1*. Click Install to continue with installation or Back if you want to change settings.

Setup - Metrel Medical ES Manager 1.22 Alpha —	×
License Agreement	
Please read the following important information before continuing.	
Please read the following License Agreement. You must accept the terms of this agreement before continuing with the installation.	
End User Licence	
Agreement	
Important Notice	
○ I accept the agreement	
• I do not accept the agreement	
Next > Canc	el
Setup - Metrel Medical ES Manager 1.22 Alpha —	×
Select Additional Tasks	
Which additional tasks should be performed?	
Select the additional tasks you would like Setup to perform while installing Metrel Medical ES Manager, then dick Next. Additional shortcuts: Create a desktop shortcut File extensions: Associate Metrel Medical ES Manager with the .padfx file extension < Back Next > Cancel Setup - Metrel Medical ES Manager 1.22 Alpha —	el
Ready to Install Setup is now ready to begin installing Metrel Medical ES Manager on your computer.	
Click Install to continue with the installation, or dick Back if you want to review or change any settings.	
Additional tasks: Additional shortcuts: Create a desktop shortcut File extensions: Associate Metrel Medical ES Manager with the .padfx file extension	
4	
< Back Install Can	cel

Figure 4.1: Survey of Metrel Medical ES Manager setup options

During the installation the Device driver Installer wizard is opened to guide you through measuring instrument driver installation. Select Next to proceed. Survey of installed drivers is presented on the screen, see Figure 4.2. Select Finish to switch back to Metrel Medical ES Manager Installation wizard. Option 'Launch Metrel Medical ES Manager' is checked by default. Select Finish to exit installation and start Metrel Medical ES Manager.

Completing the Device Driver Installation Wizard					
The drivers were successfully ins You can now connect your devi came with instructions, please re	ce to this computer. If your device				
Driver Name CDC Driver (usbser) Port	Status Ready to use Ready to use				
✓ MeasurementInstrument	Ready to use				

Figure 4.2: Survey of Measurement instrument driver installation

4.3 Metrel Medical ES Manager first use

When Metrel Medical ES Manager is first launched, Instrument pairing wizard (*Figure 4.3*) appears on the screen. It will guide you through automatic synchronization of licence purchased with Metrel instrument. Connect Metrel instrument to the PC and select Pair my instrument. Metrel Medical ES Manager will automatically collect instrument data.

Automatic pairing can be skipped on the first use and done manually at any later time when the instrument is connected to the PC.

During Instrument pairing, Metrel Medical ES Manager can find firmware upgrade for the connected instrument, provided the PC is connected to the internet.

Pairing instrument	
	Welcome to the instrument pairing wizard
	This wizard will guide you through tasks to help you to pair your instrument with the software.
	To continue, click Next
	< Back Next >

Figure 4.3: Instrument pairing wizard

4.4 Software updates

Internet connection is required for automatic detection of the latest Metrel Medical ES Manager update. When the new version of software is detected, Metrel Medical ES Manager will display a message with a link in the bottom right corner of the screen during every start-up. Click on it to download and install the latest version of Metrel Medical ES Manager software.



Figure 4.4: Update message and a link

Alternative access to software update is through the About screen. For details regarding the new software version (new features, bugs fixed, ...) see Release notes.



Figure 4.5: About screen

5 Getting started

5.1 Quick guides and tutorials

5.1.1 Introductory Guides

Before starting to use Metrel Medical ES Manager, here is some information about its basic features. The following sections will explain the way data is organized in Metrel Medical ES Manager, show you how to use advanced features and perform your work effectively.

- Metrel Medical ES Manager user interface Components
- Creating Data Structures
- Communication with the instruments.
- Measurements
- Measurement organiser
- Printing Results
- Creating Reports
- Troubleshooting

6 User Interface Components

6.1 Introduction

After start-up, Metrel Medical ES Manager Welcome screen appears with Menu tabs (●) on the top and Home tab active in the working area (②). When connected to the internet, update status presented in bottom right corner (⑤) is automatically checked.

He He	• Databa	se View	E		Home - Metrel Medical	ES Manager		
me Ner	w Open	Get Data	Connect Commun	Get instrument info nication	Measurements Organ	izer Upcomming retests Tools	Cloud	
iome ×	[4	Ð			
Dem	o files				U			
(ic\Docume	nts/Metrel	adfx Medical MESM\Dem	oFiles			
	Electrical su			Medical MESM\Dem	oFiles			
1	ntensive ca	re unit.pa ic\Docume	dfx nts\Metrel\	Medical MESM\Dem	oFiles			
(
F	Patient moni C:\Users\Publ	itor period lic\Docume	dic test de	emo.padfx Medical MESM\Dem	oFiles			
F	Patient moni C:\Users\Publ	ic\Docume	dic test de nts/Metrel/	emo.padfx Medical MESM\Dem	oFiles			
Rece	C:\Users\Publent	files	nts\Metrel\	emo.padfx Medical MESM\Dem renodic tests/intensr				

Figure 6.1: Metrel Medical ES Manager startup screen

User options available from welcome screen Home tab working area are:

- 📄 🖹 : Open existing or start new data structure file
- **Demo Files**: Open embedded demo data structure files from the list
- Recent opened files: Open data file from the list

6.2 Menu Tabs

Menus are organized in tab style, see **Figure 6.2**. Each menu tab opens a subset of functions:

- **Home** tab is active by default. It provides three groups of commands: Document for managing of data files, Communication with instruments, and Tools for starting Measurement organiser.
- **Structure** tab provides Edit commands, Structure elements and Measurements to build custom element structure test database. Testing commands enable to Run tests on the connected instrument, and View of test content.
- **Database** tab provides the client contacts data organization and the structure elements' name customization.
- View tab provides commands for main working window appearance organization.


• Main menu (the leftmost) tab provides Metrel Medical ES Manager Version data and Help file. Functions include data Import, Export and merging of data files. Settings should be checked before first start.



Figure 6.2: Menu tabs

6.2.1 Main tab menu functions

When Main tab is selected, drop down menu for function selection appears on the screen. Document handling functions are common with the Home tab menu, other functions are specific and can be accessed from Main tab menu only.





Main tab menu specific functions:




Example: Combines two Metrel Medical ES Manager Data structure files into single file. From Merge files menu, user can select master file and target file. All data are merged, same items contained in target file are overridden with master file items.



_____: Support for importing data;

Options:

Import from Exel: Import data from Excel table



6.2.1.1 Settings Menu

Settings Menu functions are organized in tab style, see Figure 6.4: Settings Tabs example.

eneral	Language	Report settings	Cloud					
)ocumer	nt path							
		Documents\Metrel	Medical N	/IESM	Brow	se		
Print p	age numbers	on report						
Print p	age numbers	on report						
Print p	age numbers	on report						
Print p	age numbers	on report						
Print p	age numbers	on report						
Print p	age numbers	on report						

Figure 6.4: Settings Tabs example

General tab functions

Set default folder for Metrel Medical ES Manager Documents. Use browser to select existing or set new folder on PC, then confirm selection.

Default document path is offered every time the document handling functions (Save, Save as, Open) are activated. It is also default search folder when Upcoming retests tool from home tab is activated. **Options:** Print page numbers on reports. Uncheck if manual page numbering is required for reports.

Language tab

• English (United Kingdom) : Current language is marked with blue dot in selection field. Select user language from the list. New selection will take effect after Metrel Medical ES Manager restart.

Report settings tab

User can choose a custom logo to appear on the printout of the results.

Cloud tab

For communication with Metrel Cloud, User API key should be generated in Metrel Cloud and copied into the API key field.

API key	*****	
Try t	o establish a communication	

Option:

Use cloud as default save location : Check to select Metrel Cloud as default save location.

6.2.2 Home tab functions

Function groups provided from Home tab menu are:

- Document
- Communication
- Reports
- Tools

a 🗈	•	B •					Intensive care unit.pa	adfx - Metrel Med	lical ES Manager		-	×
III	Home	Structu	re Datab	ase View								0
1		-	- 8	♦			0			Print Results		
Home	New	Open	Save	Get Data	Send Data	Connect	Get instrument info	Web Reports	Measurements Organizer	С	loud	
	Do	cument			Con	nmunicatio	n	Reports	Tools			

Figure 6.5: Home tab menu

6.2.2.1 Document





New : Creates a new, empty Metrel Medical ES Manager document/database.



Open : Opens an existing Metrel Medical ES Manager file from PC or Cloud.



6.2.2.2 Communication

Communication menu contains commands related to downloading data from / uploading data to instruments.



Get Data : Gets data from instrument.



Send Data : Sends data to instrument.



Connect : Opens menu window for managing transfer of group of data files between connected instrument and PC



• Get instrument info : Basic information and Firmware update status (if PC is connected to the internet) of the instrument are presented on the screen.

6.2.2.3 Reports



• Web Reports: Opens window with list of Web Report templates from Cloud. Select template and *Create report* command to open Web Report Editor application.

6.2.2.4 Tools

• Eint Results : Opens Results report in Preview window, from which page layout could be edited before printing or exporting to pdf or image file.



• Measurements Organizer : Application for creating Groups of Auto Sequence[®] tests and Custom inspections. Refer to instrument Instruction Manual for more information on Single tests and Auto Sequence[®] tests creation.



Upcomming retests
 Opens menu window for managing Scheduled retests.
 Note: Document path should be set first (see chapter 6.2.1.1 Settings Menu); Scheduled retest
 editor will automatically search for appliances within Data structure files saved in the set folder
 location.



Cloud: Opens Metrel Cloud location, if User is already signed in. If not, Sign-in window is opened.

6.2.3 Structure tab functions

Structure tab provides Structure elements, Measurements and necessary editing commands to create custom tree-like structure of the test project. Structure tab Menu groups (**Figure 6.6**) are:

- Edit
- Structure
- Measurements
- Testing
- Filtering

Elements of the tree structure can be equipped with measurements to be performed on test project. When data communication is active between PC and instrument, measurements can be performed directly from the PC on the instrument with connected test object.

Test project data file can be stored on PC and uploaded to the instrument also.



Figure 6.6: Structure tab menu groups

6.2.3.1 Edit menu group

Basic editing command icons are provided in Edit menu group. By clicking on icon, command is executed on selected data structure element. By mouse right-clicking on data structure element, context sensitive menu appears on the screen, from which additional editing commands can be selected (**Figure 6.7.**). Depending on the node selected, different menu items which apply to that type of node will be displayed.



: Copies selected data structure element together with all subtree data structure elements and associated measurements to the clipboard.



• Paste as new : Paste data from the clipboard to the selected data structure element within existing or another data structure file. The results of measurement are not included.

Measurements are pasted as new empty measurements. Command is disabled if no data was copied or if the current selected element does not allow pasting.



Paste as same: Paste data from the clipboard to the selected data structure element within another data structure file. When within same file, command is executed as Paste as new; duplicates of unique elements are not allowed within same file. Command is disabled if no data was copied or if the currently selected element does not allow pasting.



Delete : Deletes the currently selected data structure element together with all subtree data structure elements and associated measurements.



• Cut the currently selected data structure element together with all subtree data structure elements and associated measurements to the clipboard. Only active in the Results area.



• Paste cut item(s): Paste cut items from the clipboard to the selected data structure element within existing or another data structure file. Command is disabled if no data was cat or if the currently selected element does not allow pasting.

√■ • Room 01	+	
~ 🌆 PM01001	Copy item	Ctrl+Shift+C
> ഈ AP-BF_Vis_lsc ~ ₩ PM01002	Paste item as new Paste item as same	Ctrl+Shift+V Ctrl+Shift+W
> # AP-CF_Vis_PE ~ = • Room 02	Delete item	Shift+Delete
> PM02001	Rename Add/Edit Comment	Ctrl+R Ctrl+E
> Euro PM02002	Add Attachment	Ctrl+T
> Mo PM03001	Print TreeView Print All Nodes	

Figure 6.7: Context sensitive edit commands

6.2.3.2 Structure menu group

Menu group Structure contains structure elements for design of custom specific Data structure for each test project (**Figure 6.8**). With click on element it is inserted into the structure tree. Structure elements are in hierarchical relationship. Parent element type (currently selected element within tree structure) define which elements can be inserted as child elements. Those elements are highlighted, and the rest of elements are dimmed. By positioning the mouse pointer to the element, its name and description are presented.

<mark>></mark>	() ()	2		
	St	ructu	ire	

Figure 6.8: Structure elements

6.2.3.3 Measurements menu group

Measurements menu contains selection tools for adding Single tests, Auto Sequence[®] tests and Inspections to the structure elements.

Single tests : Opens list of available single tests, select one by clicking on it name. The window with settings for the selected test will open. Set parameters and limits, then confirm. Single test is appended to the selected element within Structure tree.



• Auto Sequence®s : Opens list of available Auto Sequence[®] tests (if Auto Sequence[®] Group file (*.atmpx) was opened before). Double-click on Auto Sequence[®] test to add it to the selected element within Structure tree.

Options: Browse for Auto Sequence[®] Group file. This menu is available by clicking on the triangle inside the button. First opens browser to navigate to the Auto Sequence[®] Group file default storage location; second option is to browse to alternative location.

Select Auto Sequence[®] Group file and confirm selection to open it. New Auto Sequence[®] Group file will replace existing.

Note: Only one Auto Sequence® Group file can be opened at the same time.



Inspection: Opens list of standard inspections. Select one by clicking on it name. The window with settings will open. Confirm selection. Inspection is added to the selected element within Structure tree.

Note: Custom created Inspections can be created and used as measurements within Auto Sequence[®] test, refer to Chapter *10 Measurement Organizer* for details. Note that testing by standard IEC 60601 requires the potential inspections to be custom created.

6.2.3.4 **Testing menu group**

Testing commands are active when a measurement in the Data structure tree is selected. Measurement must be attached to a corresponding structure element.



• View : Opens selected measurement in separate window for detailed overview of test flow, properties, settings, limits and detailed results if available.



Run: Opens selected measurement in separate window to run it on the connected instrument.

• Testing menu commands are also available within Data structure tree.



6.2.3.5 Filtering menu



• Filters: Opens Filter tab at the top of Results view area of the User workspace. Multiple Filter tabs can be opened simultaneously. See chapter *6.3.3.1 Results view filters* for more details.

6.2.4 Database Tab

Frequently used data can be saved to the organized tables provided in Database tab. This data can be appended to the Structure elements or to the Print results header (see Settings / Printout headers). Appearance of the Database tab menu is presented on **Figure 6.9**. Its main areas are:

- Contacts
- Structure names
- Custom Lists

⊠ I∎∙ Home Struct	ure Dat	Hos abase View	pital Equipmen	t Safety AStes	t_exe1.padfx -	Metrel ES Ma	anager Beta			- □ >
Clear Delete row Edit										
Contacts	^	Organizatio	n Name	Address	Telephone	Mobile	Fax	Email	Location n	Postcode
Client Adress of location		 Human he. 	. Biomedical	Sunny str	+386999099	+38699123		info@bio	1000	SI1000
Names Structure Names	^									
Custom Lists Appliance ID Name Group Location (Room)	^									

Figure 6.9: Database tab - Client list in Contacts group is selected

6.2.4.1 Contacts data lists

Organized tables are provided for custom Contacts data entries (Organization, Name, Address, Telephone, email etc.), which can be added to the structure elements with simple selection from drop-down list:

- Client (element Client)
- Address of location (element Location)

Edit menu options:



Clear : Clear all data entry from selected Contacts list.




Delete row : Delete selected row from the Contacts list

6.2.4.2 Structure Names data list

Structure names is a fixed list of elements which can be used to build custom tree-like structure of the test project (**Figure 6.10**). User can customize the name of the element under which it is presented within tree structure.

Structure Names data columns are:

- Picture: graphic presentation of the element
- Default name: Metrel Medical ES Manager Default Element name
- Custom name: User defined name of the element; select current Custom name and edit it.





Edit options:



Clear : Structure Names list is set to default state. After Warning confirmation, all Custom names are reset to default names.

Note:

Custom name is always used to present Element within tree structure when added to the location. By default, Custom name is the same as Default name of the element.

6.2.4.3 Custom Lists


Additional Custom data lists (**Figure 6.11**) are provided for entry of frequently used data of test objects, which could be appended to the structure elements presenting electrical equipment, with simple selection from drop-down list:

- Appliance ID (empty list by default)
- Name: name of electrical equipment (partially completed list by default, user can add new names)
- Group: name of equipment group (partially completed list by default, user can add new groups)
- Location (Room): name of the space location (empty list by default)

☐ Home Structure	Data	Hospital Equipment Safety AStest_exe1.padfx - Metrel ES Manage tabase View
Clear Delete row Edit	2	
Contacts	8	Name
o on la olo		CD / DVD player
Client		Digital TV box
Adress of location		Fax machine
		Freezer
Names	0	Fridge / cooler
Structure Names		Games console
Officiale Names	1	Gardening appliance
Custom Lists	8	Hairdryer
Subtorn Lioto		Hand dryer
Appliance ID		HiFi / Radio
Name		IEC lead
Group		Iron
Location (Room)		Kettle

Figure 6.11: Custom list Name

Edit menu options:





• Delete row : Delete selected row from the list.

6.2.5 View tab

View tab provides commands (**Figure 6.12**) to control appearance of the User workspace, where detailed content of opened data structure files is presented. By default, opened Data structure files are presented in cascaded tabs. Options are:

- Tile Vertically: opened files are tiled vertically within user workspace.
- Tile Horizontally: opened files are tiled horizontally within user workspace.
- Close all documents: close all opened data structure files including the Home tab of the user workspace.
- Reset window layout: switch to default cascaded tab presentation of opened files.
- Switch windows: opens drop-down menu for selection of active data structure file.




Figure 6.12: View tab options

6.3 Metrel Medical ES Manager User workspace

User workspace is presented on the **Figure 6.13** below. Multiple Data structure files can be opened at the same time. Their names are listed within document tab strip (\mathbf{O}), currently active data structure file tab is highlighted.

Data structure file content is presented in three main areas within User workspace, each providing special insight into the file data content, for easier navigation and data survey:

- Tree View window (❷)
- Result view area (6),

.

• Properties window (4)




Figure 6.13: User workspace

6.3.1 User workspace management

Result view area is always in the centre of user workspace and acts as container (parent) window to which Tree View window and Properties window are appended. When data structure file is opened or downloaded from the instrument, Tree View window and Properties window will be docked on the left and right side. User can manage appearance of workspace with:

- Resize Tree view and Properties window width: position mouse marker on window edge. When an arrow (<==>) appears, left click and drag it to the left or right.
- **Resize Result view columns**: position mouse marker on column on the left side edge in title bar. When an arrow (<>>>) appears, left click and drag it to the left or right.

Note:

Result view area is always expanded between Tree view and Properties window.

6.3.2 Tree view options

The tree view (see **Figure 6.14**) allows you to navigate through the data structure and displays the status of the measurements performed on elements by adding measurement status symbols to element icons.

When measurement is selected, Testing commands are offered for detailed results survey or execution on the connected instrument.




Figure 6.14: Tree View

Use buttons under the window title $(\mathbf{0})$ to navigate:

- Expand all nodes of the data structure.
- Scollapse all nodes of the data structure.
- Enter text to search... P: Find panel to search for the text string entered to the field. All string appearances are highlighted within data structure tree.

Use triangle mark in front of the element to expand / collapse subtree nodes:

- EXECUTE: Structure element subtree is collapsed: click triangle mark to expand node.
- Structure element subtree is expanded: click triangle mark to collapse node.

Right click on the structure element opens menu window (2) with commands:

- Copy, Paste as new, Paste as same, Delete, and Rename items are functions for manipulating data structure elements and measurements.
- Add / Edit Comment: Write short comment which appears as a string below selected element (3).
- Add Attachment: Add attachment to the selected element. Notification (¹) appears in the attachment column of the data structure tree. Attachments can be any type of data (picture, text file, element, etc), and it can be opened, edited or removed. Click on the notification to present a list of attached elements. Attachments are not available for measurements. Note:

Attachments are only stored in the data structure file and will not be sent to the instrument. In case the data structure file with attachments is sent to the instrument and then

downloaded into a new data structure file, attachments will be lost. Download the measurements into the same data structure file to keep using them.

- Print Tree View: opens the current Tree view (only selected nodes expanded) in Preview manager window. Preview contains all standard functions for editing printable files. It is possible to edit Header and Footer or set Background colour. It is also possible to export edited document in .pdf format.
- Print All Nodes: expands all nodes and opens the Tree view in Preview manager window.
- Instrument data: opens a window with the data of the instrument with which measurement was performed. Available for finished measurements (single tests or Auto Sequence[®] tests) only.

Testing commands buttons (④) appear at the right end of selected measurement, used for:



• View : Opens selected measurement in separate window for detailed survey of test flow, properties, settings, limits, and detailed results (if measurement is already executed).

• Run : Opens selected measurement in separate window with commands provided to set communication port and run measurement on the connected instrument.

Structure element test status indication

Auto Sequence[®] and other structure elements' test status is indicated with a dot following the icon:

- Empty status (not yet performed) measurement is associated with element; click the triangle in front to expand branch and view content of the measurement.
- Performed measurement with Fail status is associated with the element; click the triangle in front to expand node and view individual test status.
- Performed measurement with Pass status has no special indication; click the triangle in front to expand node and view individual test status.
- Element without a triangle mark in front has no measurements associated.
- The element subtree has a Fail status.
- There are measurements with Empty status and with Fail status within element subtree.

Single test measurement status indication

Measurement status is indicated with circle in front of Single test name.

- Passed finished measurement with test results performed by the instrument.
- Failed finished measurement with test results performed by the instrument.
- Finished measurement with test results and no status performed by the instrument.
- **O**: Empty measurement without test results.

6.3.3 Results view area options

Table with extended information organized in columns is located within the Result view area (Figure 6.15):

- Appliance ID: identification of the tested appliance.
- Structure path: Appliance directory within the data structure tree.
- Name: Appliance name.
- Location (Room): Appliance location
- Next test: Date of next test for the Appliance with set Test period.
- **Test date:** Test Date from the instrument.
- Status: Test status

The Appliances table is focused on the selected element within the Tree View; measurement results associated with the appliances located within the subtree structure of the selected element are listed in the table.

Interactively, when a left click is applied on the Appliances table row within the Result view area, an adequate Tree view branch is expanded, and the selected item is highlighted.

Appliance test status and appended **measurement status** is written in Status column; click on + in front of table row to expand content and survey measurement status details:

- **Pass**: Passed finished measurement with test results.
- **Fail**: Failed finished measurement with test results
- Nothing: Finished measurement with test results and no status.
- **Empty**: Empty measurement without test results.

Applian	nce Id 🔺	Structure Path	Name	Location (Room)	Next test	Test date	Status	
PM0	1001	Demo PM Saf	SenseVis BF50	IC Room 1	20/04/2021	20/10/2020	Pass	
ΞA	AP-BF_Vis	s_lso_Leak					Pass	
	Visual	IEC/EN 62353					Pass	
	Riso L	N-NEP					Pass	
	Riso L	N-NEP					Pass	
		Pass						
	Pass							
	Riso L	N-NEP					Pass	
	Applied	d Part Leakage (dir	ect)				Pass	
PM0	1002	Demo PM Saf	SenseVis CF80	IC Room 1	20/04/2021	20/10/2020	Pass	
B AP-CF Vis PE Leak								
		Pass						
	Pass							
Protective Earth Resistance								
	Protec	tive Earth Resistan	ice				Pass	
	Applied	d Part Leakage (dir	ect)				Pass	
PM0	2001	Demo PM Saf	SenseVis BF5	IC Room 1	20/04/2021	20/10/2020	Empty	
ΞA	AP-BF_Vi	s_lso_Leak					Empty	
[Visual	IEC/EN 62353					Empty	
	Riso L	N-NEP					Empty	
	Applied	d Part Leakage (dir	rect)				Empty	
PMO	2002	Demo PM Saf	SenseVis CF80	IC Room 1	20/04/2021	20/10/2020	Empty	
PM0	3001	Demo PM Saf	SenseVis BF5	Surgery	20/04/2021	20/10/2020	Empty	
PMO	3002	Demo PM Saf	SenseVis CF80	Surgery	20/04/2021	20/10/2020	Empty	
gennerrer	DOF M	s PE Leak					Empty	

Figure 6.15: Results view area

Note:

Appliance table within Results view area present status of listed Appliances, Auto Sequence[®] tests, Single tests and Inspections. For detailed survey of test settings, results, sub-results, parameters, limits use Testing command View.

6.3.3.1 Results view filters

Sort filter

• Test date •: Click on the column name to activate sort filter indicated with arrow.

Each Appliances table column provides a sort-by-name filter. Active sort column is marked with arrow on the left side of the column name field. Click again to reverse sort order. Sort filter can be active in only one column at the same time.

Appliance table column filter

• Name : Position mouse over column name field and click filter sign appearing on the right-up corner of the cell to activate column filter.

Each column of the Appliance table contains a filter. After filter is selected, menu with filter options appears on the screen. Select the option and close the window to activate the filter. Active column filter remains marked with the blue filter sign within name field, see **Figure 6.16** below.

Appliance Id 🔺	Structure Path	Name 📍	Location (R	Next test	Test date	Status	Ŷ
E PM02002	Demo PM S	SenseVis CF	IC Room 1	20/04/2021	20/10/2020	Empty	
# AP-CF_	Vis_PE_Leak					Empty	
E PM03002	Demo PM S	SenseVis CF	Surgery	20/04/2021	20/10/2020	Empty	
AP-CF	Vis_PE_Leak					Empty	

Figure 6.16: Results view area column filters – Name and Status are active.

More column filters can be active at the same time.

To remove column filter, select it and choose the Clear Filter command from filter settings menu.

Appliances table Custom Data filter

Custom Data filter can be active within any column of the Appliances table, with one or more filters active at the same time.

- Custom filter: Right click within Results view area to open a menu window and check Custom filter option to open Custom filter design tab at the top of Result view area (•), see Figure 6.17.
- Position mouse over any field within filter tab and click to open additional filter criteria or change or delete it.
- Options available for manage designed Custom Data filter are:

Apply filter: Filter is applied to the Appliance data table, active data filters are marked within data column name field.

Save filter: Menu is opened to Name and save filter in Filters folder Medical MESM application. **Load filter:** Load filter from Filters folder of Medical MESM application.

Save filter as: Menu is opened to Rename and save filter.

New filter: Opens new filter tab with logic function relationship field by default.

Row filter
 Custom filter design is offered as a row beneath Address row of the Appliance table
 (2), see Figure 6.17. Custom filter and Row filter are equivalent. Settings of the Custom filter are mirrored to the Row filter and vice versa.

	Status	nce Id Begins = Empty est < 08/05/20	1							
	Apply	filter	Load filter	Save filte	ər	Save filt	er as	New	/ filter	
A	polian 🔺 📍	Structure Path	Name	Location (R	Next	test 📍	Test date	2	Status	
-	R PM	=	= 2	REC .	0.000000	/05 🗸	=	~	= Empty	
Ξ	PM02001	Demo PM S	SenseVis BF	IC Room 1	20/04		20/10/20	20	Empty	
	AP-BF_V	/is_lso_Leak							Empty	
and the second se		SenseVis CF IC Room 1		20/04/2021 20/10/2		20/10/20	20	Empty		
	H AP-CF_V	/is_PE_Leak							Empty	
Ξ	PM03001	Demo PM S	SenseVis BF	Surgery	20/04	/2021	20/10/20	20	Empty	
	AP-BF_V	/is_lso_Leak							Empty	
Ξ	PM03002	Demo PM S	SenseVis CF	Surgery	20/04	/2021	20/10/20	20	Empty	
	AP-CE V	is PE Leak							Empty	

Figure 6.17: Appliance Custom data filter design options

Appliances table Find filter

- Find panel: Right click within Results view area to open a menu window and check Find panel option. Field for Find filter is opened at the top of Result view area (**0**), see Figure 6.18.
- Filter the table to only rows that contain the chosen string (2). Search is started when Find button is clicked.
- Appliances table is filtered to the rows containing the string. (3). The string appearances are highlighted.
- Since the second seco

BF50	9	8 × 1	Find		
Appliance ID	Structure Path	Name	Next test	Test date	Status
E PM01001	Demo PM Safety testing/General Health	SenseVis BF50	20/04/2021	20/10/2020	Pass
⊞ AP-BF_V	ïs_lso_Leak				Pass
	Demo PM Safety testing/General Health	SenseVis BF50	20/04/2021	20/10/2020	Empty
E PM03001	Demo PM Safety testing/General Health	SenseVis BF50	20/04/2021	20/10/2020	Empty
AP-BF V	is Iso Leak	0			Empty

Figure 6.18: Appliance data Find filter design options

Appliances table commands

Right click within Appliances table offers a menu window from which commands can be selected:

- Clear filters : Clear all filters from Appliances table.
- Delete item : Delete selected item.
- Delete all : Delete all items from appliances table.

6.3.4 Properties window options

Properties window contains description of any element within the project data structure tree. Properties window appearance depends on the properties of the active element type.

The Properties window view in Figure 6.19 is focused on an element selected in the Tree view.

- Selected element properties and properties of all parent elements within the data structure are presented in separate tabs in the Properties tab strip (1). By default, selected element tab is active. Click on the parent element tab to see its properties. Use arrows at the end of the tab strip to scroll between tabs.
- List of available properties to the element is presented in the left column (♥). Click the entry field on the right side and to input the property value.
- When a drop-down list of predefined properties is available, arrow will appear at the end of entry field (3). Click on the arrow to open drop-down menu and select property value.




Figure 6.19: Properties window

7 Creating Data Structures

7.1 Considerations on a new test project's data structure

Metrel Medical ES Manager supports creating a tree-like data structure of a test project, with objects under test located on the tree nodes in a hierarchical relationship. Measurements can be appended to the objects before uploading the data structure file to the instrument.

Before opening a new data structure file of a test project, consider:

- **Data structure tree**: determine the data structure of the test project test objects and their hierarchical relationship within tree-like project presentation.
- Data lists with properties: consider the test elements properties needed for consistent results presentation and the test report. Commonly used data can be saved to the lists in the Database tab.
- Auto Sequence[•] tests: For faster test execution, collect the tests necessary for each individual element and create an Auto Sequence[®] Group of tests. Auto Sequence[®] test can contain custom created Inspection checklists. Refer to the individual Instrument Instruction manual for details.

7.2 Creating a new data structure file for a test project

To create a new data structure file for a new test project, follow the procedure:

- Optionally, set the document path from Main tab \rightarrow Settings \rightarrow General.
- Select Home tab of the user workspace.
- Chose New file from Home user workspace or Home tab (or Main tab).

New data structure file *.padfx is opened, with default file name *New Document**. User Workspace area is empty, except the parent element *Node* is automatically placed into the new Tree view. This is the starting point to create the project data structure. It cannot be deleted: one *Node* must always exist in the structure tree. It can be renamed to reflect the test project subject.

Before starting to create a data structure it is recommended to:

- Rename the file using the Save as command from the Home tab menu, to reflect the test project subject.
- Open the appropriate Auto Sequence[®] group file containing the required project tests, if available. Only one Auto Sequence[®] group can be opened at the same time.

7.2.1 Creating a new data structure

Creation of a new Data structure starts from the tree view of the user workspace, see **Figure 7.1**.

- Select Structure tab menu (**0**). Available elements are highlighted within Structure group menu.
- Default parent element Node is already selected within the Tree view. Click on the element in the Structure menu group to insert it into the structure tree. It appears as a new tree element (2), hierarchically as a child element under the parent Node. To place another tree element, repeat the action.
- To add more elements and levels to the tree structure, select an element in the tree view and click on a new element from Structure menu group. It appears as a child element under the selected one. To further expand the levels, select the child element and repeat the action.

- To create another node subtree, select the parent node and repeat placing new elements as child items on the same level (③). The parent element and its children form a subtree.
- Describe the elements added to the tree structure in the Properties window of user workspace

 (④). When Name property is entered, it is automatically updated in the Tree view as well. Each element has a dedicated set of properties, which can be manually entered or selected from the embedded or user created Database lists. Describe all the properties necessary for the desired level of detail for interpreting test results and for well-defined Test Reports.
- Repeat the above actions until the desired test project data structure is finished and save the file.

	New D	ocument.padfx - Metrel Medical ES Manager		×
Home Structure Database View Structure Database View Copy Paste as new Paste as same Delete Cut Pated Edit	Structure	Single tests Auto Sequence®s Inspection Measurements	View Run Filters Testing Filtering	Ø
Home × Hospital Equipment Safety AStest.padfx ×	New Document padfx ×			4
Tree View	Appliance Id Structure Path Nar	ne Location (Room) Next test Tes	st date Status	Properties
Enter text to search	Medical e Node/Locati		Empty	S Location1
Name	Medical e Node/Locati Medical e Node/Locati		Empty	General
	Medical e Node/Locati		Empty	Name (designation) of Location1
✓ ≥ Node ✓ SLocation1	-			Description of location
✓ Project1 ✓ Element1				✓ Adress of location ✓ Organization Name Address Telephone Mobile Fax Email Location number Postode
 ♥ Project1 ♥ Element1 ♥ Medical equip.1 ♥ Element2 ♥ Element2 ♥ Medical equip.1 ♣ Medical equip.1 	Row count: 4			Postcode

Figure 7.1: Creating a new Data structure

7.2.2 Appending measurements to the data structure elements

It is possible to prepare required measurements (Single tests, Auto Sequence[®] tests and Inspections) and append them to the test elements.

- The Data structure with pre-prepared measurements can be then uploaded to the instrument.
- Another option is to start measurement on the instrument directly from PC's Data structure.

7.2.2.1 Appending Single tests



- Click on icon Single tests opens a list of available Single tests, organized in groups.
- Expand single test groups and select the single test from the list by clicking on it name. Single test setting window appear on the screen, see chapter *9.1.2 Viewing or configuring the single* tests for setting detailed information.

7.2.2.2 Appending Auto Sequence[®] test

To append an Auto Sequence[®] test, the Auto Sequence[®] Group file that contains it should be opened first. Auto Sequence[®]s command from the Measurements group of Structure tab menu has two options:

Browse for Auto Sequence® group file:



- Click on the triangle Auto Sequence®s opens a menu with the command for browsing for the Auto Sequence® Group file. Select the file and confirm selection to open it. The new Auto Sequence® Group file will replace the existing one. Only one can be opened at the same time.
- Appending the Auto Sequence[®] test:



- Click on the icon Auto Sequence®s opens a list of available Auto Sequence® tests from active Auto Sequence® Group file.
- Expand the folders to find the requested Auto Sequence® test.
- Double-click on the Auto Sequence[®] test to automatically append it to the currently selected element.
- The set connections, parameters and limits of the Auto Sequence[®] test can be reviewed. Select Auto Sequence[®] test and click on the icon ^Q on the right side of the tree view to open survey window. See chapter *9 Measurements* for more information.



Procedure: Append Inspection test Inspection

- Click on the icon opens a list of available standardized Inspections, organized in groups.
- Expand the Inspection group and select the Inspection from the list by clicking on it name.
- Inspection window opens with the default Empty status of inspection items.
- At the top of the tab the meta-data (date and time) can be added to the configurable inspections.
- Click OK to confirm. Inspection is appended to the selected element.

Note:

Only the standardised Inspections are available from Inspection menu list.

Custom Inspections can be created and used as measurements within the Auto Sequence[®] tests. Refer to Chapter *10 Measurement Organizer* for details.

7.3 Creating Data structure using Copy and Paste commands

If parts of the data structure are often repeated, it is possible to create large structures with measurements in a fast and simple way with the help of Copy and Paste commands.

Objects from one data structure can be copied to the structure in another file. Data structures can be copied without measurements, with appended empty measurements or with executed measurements.

Note

Each element and measurement inserted to the data structure has got an internal unique identification. Only one item with same identification can exist within same data structure. See below for more information.

Commands are available from the right-click menu on the element in the Tree view, or from Structure tab Edit menu:

• Copy command:

- Selected element and all elements located within the selected subtree structure are copied.
- o Measurements appended to the element, empty or executed, are copied.
- Properties of each copied element are copied.
- Paste as new:
 - Copied element and all element located within it subtree are pasted to the new selected location.
 - Empty measurements appended to the elements are pasted and remains unchanged.
 - o Executed measurements are pasted as empty measurements.
 - All pasted elements and Measurements got new internal unique identification.
 - Properties of each element are pasted and remain unchanged.
- Paste as same:
 - Copied element and all elements located within its subtree are pasted to the new selected location.
 - Measurements appended to the elements, empty or executed, remain unchanged.
 - o All pasted elements and measurements keep the same internal unique identification.
 - o Properties of each element are pasted and remain unchanged.
 - o Only possible when pasting into a new file.

Note:

When Paste commands are dimmed, pasting data to the selected location is not allowed due to hierarchical relationship between the copied element and paste element, or no data was copied before.

7.3.1 Example of Copy and Paste within the same data structure file

Copy and Paste as new: a structure element with a sub-tree and either with empty measurements or without them.

Command sequence is useful to expand the data structure with the repeated common nodes, like multiple units, each containing equally equipped rooms (Figure 7.2):

- Select the node with a subtree, for example Intensive care unit (**0**) and execute Copy command.
- Select the new location, for example its parent node Patient Monitor Safety testing (2) and execute Paste as new command; repeat for each new node.
 - Paste as new was repeated twice (③), creating two new nodes Intensive care units.
 - All elements within pasted nodes have new unique identifications, but it names and Properties and appended measurements remains the same.
- Select each new element within new nodes and update only their properties (④), which distinguish them from the copied node elements.




Figure 7.2: Copy and Paste as new operation within same Data structure - empty measurements appended

Copy and Paste as new: a structure element with a subtree containing finished measurements. The command sequence creates a new data structure with the same elements with the same properties, while the measurements paste as empty. It is useful to perform the same test project again on same or another equal location (Figure 7.3):

- Select the root of the subtree containing finished measurements (**1**) and execute the Copy command.
- Select a new location (the new Node in example below) (2) and execute Paste as new command:
 - Pasted data structure tree and element properties remain the same. All elements have a new internal unique identification to distinguish them from the source elements.
 - All appended measurements are empty (③), without results. Parameters and limits settings remain the same as they were in the copied measurements. All measurements have new internal unique identification.





Figure 7.3: Copy and Paste as new operation within same Data structure - finished measurements with results

Copy and Paste as same:

Internal unique identification rule prevents Paste as same command execution within same data structure file. Paste as same command is executed as Paste as new in that case. An information message (Figure 7.4) appears on the screen.





7.3.2 Copy and Paste between different data structure files

Source data structure file and the target data structure file should be opened at the same time.

Copy and Paste command sequence execution always begin with same procedure:

- Select source file tab to become active within user workspace.
- Select the root of the chosen subtree from the Tree View and execute the Copy command.

Paste as new:

- Select target data structure file tab to become active within user workspace.
- Select parent element in Tree view in the target data structure and execute the Paste as new command.
 - Pasted Data structure tree and **elements'** properties remain the same as the source. All elements have new internal unique identifications to distinguish them.
 - All appended measurements from the source file, empty or done, are pasted as empty. Parameters and limits settings remain the same as they were set in the source file. All measurements have new internal unique identification.

The difference between Copy and Paste as new between different files is only in presence of source part of the tree in the target file. For examples of use refer to the chapter **7.3.1 Example of Copy and Paste within the same** data structure file.

Paste as same: Source data structure elements have appended empty measurements or no measurements:

- Select target data structure file tab.
- Select the chosen parent element in the tree view and execute the paste as same command.
- Unique identifications of elements and measurements are automatically compared.
 - Same unique identifications were found: merge confirmation information window on Figure 7.5 appears on the screen.
 - Click OK to complete the Paste command with data merging. Another confirmation is needed.
 - o Elements and measurements keep internal unique identification.
 - o Click Cancel to stop command execution. Data will not be pasted.




Figure 7.5: Merge confirmation information window

- **No equal unique identifications**: procedure continues with the confirmation window, **Figure 7.6**.
- Click Yes to complete the paste command. Elements and measurements keep internal unique identification.
- o Click No to stop command execution, data will not be pasted.




Figure 7.6: Paste as same confirmation information window

Paste as same: source data structure elements with appended finished measurements; target data structure has empty measurements, or finished measurements executed with different instruments:

- No equal unique identifications within source and target data structures:
 - o Procedure is the same as described above.
 - Equal unique identifications within source and target data structures:
 - Procedure starts the same as described above until merge confirmation step.

Paste as same: source data structure elements only; target data structure elements only:

- User is prompted to select next paste command step, Figure 7.7:
 - **Merge duplicates using source data**: target data structure elements with equal identification are updated with source data. The rest are pasted the same as source data.
 - **Relocate structure**: source data will be pasted to Target Data structure as additional node subtree. New elements have new internal unique identification.
- Conclude paste command with response to Information message on the screen, already described above.

Download data
Download to New file
Download to Currently opened File
Cancel

Figure 7.7: Paste as same confirmation information window

7.3.3 Example of merge functionality

We have a test project with multiple identically equipped locations (rooms) and we require a test report for the complete structure.

Complete structure and empty measurements of the test project are uploaded to multiple instruments. Each tester measures its dedicated locations only. Data from all instruments is downloaded to the same Metrel Medical ES Manager.

If the test report is created from the each individual data structure files, each will include complete structure but only the elements downloaded from the relevant instrument will contain finished measurements. This system doesn't fulfil the requirements.

Merge function enables all the downloaded data structure files from different instruments to be merged into a single file that contains the complete structure and complete measurements. Complete test project report can now be created.

Special case: Source data structure elements with appended empty measurements; target data structure elements with appended finished measurements:

• Merged data structure and elements remain the same, empty measurements are appended to the appropriate structure elements in addition to the unchanged finished measurements.

7.3.4 Example of copy and paste as same functionality

We have a test project with multiple identically equipped locations (rooms) and we require a test report for the complete structure.

Same data structure that covers only each single location with measurements is uploaded to multiple instruments. Each tester measures its dedicated location only. Data from all the instruments is downloaded to the Metrel Medical ES Manager.

If the test report is created from the each individual data structure files, will include only **the each instrument's** measurements and structure, which is not what is needed.

Copy and paste as same function enables data from all locations to be combined into a single document file. Complete test report can now be created.

8 Communication with the instruments

8.1 Establishing communication with instruments

Communication with instruments is established via RS-232 port or USB port.



- Click on Get Data, Send Data or Get instrument info to open a dialog box for setting the parameters of communication.
- From drop-down list select COM port to which instrument is connected.

Serial	Ethernet				
Basic	Advanced				
COM	port				
COM	3: Measurem	ent Instrument	USB VC	1	

Figure 8.1: Setting communication

8.1.1 Instrument info

Basic information of the instrument (name, short name, manufacturer, type of instrument, serial number, HW version, FW version, calibration date, etc.) and firmware update status can be obtained. If no data is displayed, see the chapter *13 Troubleshooting*.

Options:

- Check for new updates: check if new FW version exists
 Download: Download the new instrument FW version
- View release notes: open release notes related to the instrument type

Name	MediTest
Short Name	MI 6601
Serial	21321253
Firmware version	1.2.22.3e3983fd
Instrument code	CAAA
Hardware version	1
Calibration date	06/09/2021
Manufacturer	
	1
re are no new updates ava	ilable.

Figure 8.2: Instrument info screen options

8.1.2 Upgrading the Firmware

It is possible to upgrade the firmware on the instruments to the latest version using the Metrel Medical ES Manager. Internet connection is required. When the instrument info window notifies about an upgrade, click Download and follow the instructions. Upgrading program (Flashme) will start and guide you through the upgrade procedure.



Note

For detailed set up procedure please see the user manual of the attached instrument on your computer! User manuals can be downloaded from our <u>Download centre</u> [http://www.metrel.si /download-center.html].

8.2 Upload or download Project files

Upload or download of the data structure project files are located in the communication group of the Home tab menu. Single and multiple file transfer is supported.

8.2.1 Single file transfer options



Get Data: Gets data from instrument. The currently active data structure file in the instrument memory organizer is downloaded.

- If only the Home tab of the user workspace is active, the downloaded file is opened in a new user workspace tab.
- If the data structure file is already opened in the user workspace, download options selection window is opened, see Figure 8.3:
 - **Download to a new file**: new data structure file is opened in a new tab in the user workspace.
 - **Download to the currently opened file**: elements data with equal unique identification are merged, rest of them appear as new elements within the active data structure file.

Download data	
Download to New file	
Download to Currently opened File	
Cancel	

Figure 8.3: Download options selection box



- Send Data: Upload the active data structure file to the instrument.
- 8.2.2 Multiple file transfer options

- Connect : Open the menu window (Figure 8.4) for managing the transfers of multiple data structure files between the connected instrument and the PC. The window contains Transfer tab and Instrument info tabs.
- Transfer options are:
 - o **1** Serial communication port selection and the Reload command.
 - PC folder selection tool. Use the arrows to navigate the folder structure, or open the folder content where the project data structure files are saved.
 - Elist of the data structure files in the selected folder on the PC. Click on the square in front of the file name to select the files for upload to the instrument.
 - List of the data structure files in the instrument memory. Click on the square in front of the file name to select the files for download to the PC.
 - Command options:

0

0

- Upload selected Project data structure files to the instrument.
 - Download selected Project data structure files to the PC folder.
- o Delete selected Project data files from the instrument.
- o Connection Select wired connection Serial or Ethernet.
- o Select the PC serial data port that the instrument is connected to.
- o Host name / IP Enter IP address and Port number
 - Reload
 Reload/refresh the data.

	Connection - 🥝 Reloa	d 😭 Port COM3:	Measu	rem 🚩	
Transfer	1 Instrument info				U
-	↑ ** M → Strukture in AS in cus → × → 2		D	elete sele	ected
	Name				Name
¥	Hospital Equipment Safety AStest				piezo zaga
~	Hospital Equipment Safety AStest_exe1				Workspace_001
V	Hospital Equipment Safety AStest_exe2			>	Workspace_002
	Intensive care unit				Workspace_003
	New Document				Workspace_004
	New Document_test	\bigcirc			Workspace_005
	1				Workspace_006
		0			Workspace_007
					Workspace_008
					Workspace_009
					Workspace_010
					Workspace_011
					Workspace_012

Figure 8.4: Transfer tab menu of the Connect window

• Instrument info tab: Instrument info list, see chapter 8.1.1 Instrument info.

9 Measurements

When the connection with the instrument is active, measurements can be executed from the PC using the Metrel Medical ES Manager.

9.1 Single tests

9.1.1 Single tests in the user workspace

Options:



 Single tests New single test: click on the Single tests icon in the measurement group of the Structure tab to open a list of available single tests. Select one by clicking on its name. It will open the configuration window for configuring, saving and/or executing the single test on the instrument.



• Run Run single test: Select the single test already added to the data structure tree and click the Run command found in the same line in the tree structure, or in the Testing menu tab. The test configuration window will open where the single test can be surveyed, reconfigured or executed on the instrument.



• View View single test: Select the single test from the data structure tree and click the View command found in the same line in the tree structure or in the Testing menu tab. The configuration window will open. The single test can be surveyed only, reconfiguration and execution are not possible.

9.1.2 Viewing or configuring the single tests

Single test configuration window is presented on **Figure 9.1**. Its functions are:

- **Onnections** are the test object settings:
 - o Operating parameters.
 - o Connections types.
 - o Connections functions.
- **2** Single test: parameters and limits settings.
- **B Results:** detailed list of the single test results.
- **4** Summary: summary of the results state.
- **Commands:** running the test, and the communication port survey.

© Conn	ECTIONS	Selected Port COM7	Speed 115200	- E	
PARAMETERS		RESULTS			
Class	1.	~ Patient Leakage (to earth	Ĩ.		
FE	Yes 🔻	5.0		8	
SIO	No 🔻	VMn Ffpe FEo Ne Rt L	21	9	no value
Test On delay	55 -	~ 1 VMn Ffpe FEo Ne Rt L			no value
Manual power off	Off ▼	✓ I VMn Ffpe FEo No Rt L			no value
Power Off delay	5 s v	~ I vMn Ffpe FEo No Rt L	P3		no value
CONNECTIONS		✓ I VMr Ffpe FEo Ne Rt LF	1		no value
► BF AP BF		~ 1 VMr Ffpe FEo Ne Rt LF	3		no value
NEP NEP		~ 1 VMr Ffpe FEo No Rt LF	и		no value
		~ 1 VMr Ffpe FEo No Rt L	3		no value
+ ADD CONNECTION					
AP BF	NEP				
SF1 SF2					
	00000	1 I			
1 2 3 4 5 6	5 7 8 9 10 P/S	SUMMARY		•	
	r	Elapsed time		Passes	
atient Leakage (to earth)		Started at		Fails	
Singe	E TEST	Estimated time	00:01:40	Remaining	
				Total	

Figure 9.1: Single test configuring window

Refer to the test object's instruction manual for the Connections setting details. Refer to the test instrument's user manual for the single test settings details. After the single test configuration is finished, the options are:

- Save the single test to the previously selected data structure location and close the window.
- Cancel Discard changes, the single test remains unchanged.

9.1.3 Running single tests

Open the single test configuration window and set the communication port where the instrument is connected, then start the test execution using these commands (on the Figure 9.1):

- Open the survey of communication port setting, see the picture below.
- Check the settings are correct.
- Start the test on the instrument.

During the test progress, data is acquired and presented in the Results and Summary areas in real time. Notifications and test flow commands are provided to the user, see **Figure 9.2**:

- **Results**: Single test name is presented at the top. Individual results with statuses are presented in real time, the user can scroll through results sub window during test.
- **Summary**: the test statistics are presented in real time
- OPProgress bar and percentage of the single test is calculated and presented in the top right corner of the Commands area and next to the single test name.
- **ONDER INTEGRATIONS** OF test flow appear at the bottom of Results sub window. Some of them require user action (e.g. turning the tested device on).

- **6** Commands for user test flow control:
 - Next step: continue with the Single test after the user action is completed.
 - Next test: immediately stop the Single test
 - o Stop test: immediately stop the Single test

RUNNING MEASUREMENTS 67%					
ESULTS	💲 TOTAL PATIENT LEAKAGE (TO EARTH)				
Total Patient Leakage (to earth)			67%) 0		
- I VMn Ffn FEo Ne Rt LP1			0.000 mA •		
∼ VMn Ffn FEo Ne Rac LP1		U	0.000 mA 🔹		
✓ I VMn Ffn FEo Ne Rdc LP1			0.000 mA 🔹		
✓ I VMn Ffn FEo No Rt LP1			0.000 mA 💿		
VMn Ffn FEo No Rac LP1			0.000 mA 🔹		
I VMn Ffn FEo No Rdc LP1			0.000 mA 🔹		
✓ I VMn Ffpe FEo Ne Rt LP1			no value i C		
✓ I VMn Ffpe FEo Ne Rac LP1			no value 🗢		
✓ I VMn Ffpe FEo Ne Rdc LP1	no value - C				
✓ I VMn Ffpe FEo No Rt LP1	no value. O				
		4 switch o	FF THE DEVICE UNDER TEST!		
JMMARY					
lapsed time	00:01:20	Passes	12		
itarted at	11:47:56	Fails	0		
stimated time to finish	00:00:19	Remaining	6		
		Total	18		

Figure 9.2: Single test progress data and notifications



Note

While the test is running, its window cannot be closed.

9.1.4 Single test results

Finished Single test screen (Figure 9.3) provides the following information and options:

- Overall test state and status
- Single test completion percentage and status
- BResults: single test individual results, each with status; scroll up/down to survey all.
- **4** Summary: single test statistics.
- **O**ption commands.

© CONNECTIONS		>			
PARAMETERS		RESULTS			
Class	1.	Total Patient Leakage (to earth			2 10
FF	Yes 🔻	 Total Patient Leakage (to earth)		
SIO	No 🔻	VMn En Eto Ne Rac LP1			0.000 1114
Test On delay		VMn Fn FEo Ne Rdc LP1			0.000 mA
	5 s 🔻	~ I VMn Fn FEo No Rt LP1			0.000 mA
Manual power off	Off ▼	~ 1 VMn Fn FEo No Rac LP1			0.000 mA
Power Off delay	5 s ¥	~ 1 VMn Fn FEo No Rdc LP1			0.000 mA
CONNECTIONS		~ I VMn Ffn FEo Ne Rt LP1		6	0.000 mA
BF AP BF		✓ I VMn Ffn FEo Ne Rac LP1		•	0.000 mA
NEP NEP	۲	~ VMn Ffn FEo Ne Rdc LP1			0.000 mA
+ ADD CONNECTION		~ 1 VMn Ffn FEo No Rt LP1			0.000 mA
T ADD CONNECTION		VMn Ffn FEo No Rac LP1			0.000 mA
AP BF	NEP	~ 1 VMn Ffn FEo No Rdc LP1			0.000 mA
SF1 SF2		~ VMn Ffpe FEo Ne Rt LP1			0.000 mA
	000000	VMn Efne EFo Ne Rec I D1			0.000 mA
1 2 3 4 5	6 7 8 9 10 P/S				
	•	SUMMARY			
Total Patient Leakage (to earth)		Elapsed time	00:00:58	Passes	
-		Started at	14:08:31	Fails	
Sinc	GLE TEST	Finished at	14:09:30	Remaining	
				Total	

Figure 9.3: Single test results screen

Options after the single test execution is finished:

- Save the single test to the previously selected location in the data structure and close the window.
- Cancel Discard results; user is prompted to confirm action by selecting the option from the confirmation window, **Figure 9.4**.



Figure 9.4: Discard single test results confirmation

Results after repeated single test

When existing finished single test is opened from data structure and executed again, it is saved at the same location as a new finished single test a after confirmation. The existing test remains unchanged.

Results after interruption

When Single test execution is interrupted with one of the available commands, the result screen (**Figure 9.3**) provides the information as it was at the moment of interruption:

- U The overall test state is Finished with the appropriate status information at interruption.
- Ø Single test progress percentage is less than 100 at the interruption moment.
- Results: single test individual results are presented with finished status or empty status, if not yet executed at interruption.
- O Summary: single test statistics show the number of Passed (or Failed) individual tests and number of Remaining tests with the estimated time to finish at interruption.

9.2 Auto Sequence[®] tests

9.2.1 Auto Sequence[®] tests in the user workspace

To open the Auto Sequence[®] in the user workspace it must be first appended to the data structure element. Once placed in the data structure, select it to view the options:



• Run Click on the Run command in the Testing group of the Structure tab menu or in the Auto Sequence[®] line in the tree view. It opens the Auto Sequence[®] configuring window where can be surveyed, reconfigured and executed on the instrument.



• Click on the View command in the Testing group of the Structure tab menu or in the Auto Sequence[®] line in the tree view. It opens the Auto Sequence[®] in the user workspace. It can be surveyed only, reconfiguration and execution are not possible.

9.2.2 Viewing and configuring an Auto Sequence[®]

Before execution the Auto Sequence[®] content can be surveyed. Connections and the single tests' parameters and limits, and the flow commands, can be reconfigured before execution. Auto Sequence[®] configuration window is presented in **Figure 9.5**. Its functions are:

• **Properties:** Basic Auto Sequence[®] information

Select Test optimization to organise the Auto Sequence[®] steps into the fastest sequence.

- **O** Connections: click on the field to open the test object settings menu. Setup the test object's features:
 - o Operating parameters
 - Types of connections
 - o Connection functions
 - 2

0

- 2 List of single tests and flow commands in the Auto Sequence*:
- o Click on the SINGLE TEST field to open the parameters and limits settings window.
- Click on the OPERATION AFTER TEST field to open the flow command settings window.
- Results: List of all the results in the single tests of the Auto Sequence[®].
 Click on the test name to display the detailed individual results.
- **O** Summary: summary of the results status. Individual results of all single tests are included in the count.
- **Ommands:** Run test and Communication port survey.

Refer to the test object's instruction manual to correctly set up the connections. Refer to the test instrument's user manual for the Single test details and settings.

Commands at the bottom of the Auto Sequence® configuration window are:

- Save Auto Sequence[®] with the new settings to the location in the data structure and close the window.
- Cancel : Discard the changes, Auto Sequence[®] remains unchanged.

AP-BF_Vis_PE_Leak	Selected Port COM7 Speed 115200	6	
	RESULTS		
Test optimization	 Visual IEC/EN 62353 		
Auto Sequence®code	~ Protective Earth Resistance		
B101	~ Protective Earth Resistance		
Description Safety test Class I patient monitor with BF applied parts	~ Protective Earth Resistance		
	 Applied Part Leakage (direct) 	B	
	~ Lywn Lxt	•	no value O
CONNECTIONS U	VMn LX4		no value O
Visual IEC/EN 62353	~ I VMn LX7		no value O
VISUAI IEC/EN 02303	~ I VMr LX1		no value O
OPERATION AFTER END OF TEST	✓ I VMr LX4		no value 🔿
OPERATION AFTER END OF TEST	VMr LX7		no value 🔿
Protective Earth Resistance			
Single test			
Operation after end of test	SUMMARY		
×	Elapsed time	Passes	
Applied Part Leakage (direct)	Started at	Fails	
Single test	Estimated time 00:	00:44 Remaining	
OPERATION AFTER END OF TEST		Total	

Figure 9.5: Auto Sequence® configuration window

9.2.3 Running Auto Sequences[®]

Auto Sequence[®] test configuration window must be opened with Run command and the communication port with the instrument must be set before running the Auto Sequence[®].

The progress of the Auto Sequence[®] test during the test execution can be viewed or controlled.

- Results and statuses of the tests that already finished.
- Momentary test summary.
- The appropriate commands for controlling the execution of the Auto Sequence[®] test are always available. The user is notified if any action is needed.

The test results window is active during the test (scrolling up/down to view all results is possible).

Any individual Auto Sequence[®] test step can be removed right before starting the test.

- X Remove the step with the sign at the upper right corner.
- Undo Step removed.
 Undo prompt appears on the screen for short time, click on it to cancel step removal.

The parameters in the tests and the test object connections can be edited. All changes can be saved using the Save icon that appears at the first change, if Auto Sequence[®] Group file is opened.

Save icon in upper left corner, click to save changes to currently opened Auto Sequence[®]
 Group file.

Start the execution of a single test by using the commands (on Figure 9.5):

- 2 Open the survey of communication port setting, see picture below.
- Check the settings are correct.
- Start the test running on the instrument.

The data is acquired and displayed in the Results and Summary areas in real time during the testing. Notifications and relevant test flow commands are provided to the user, see **Figure 9.6**:

• Results: List of single tests contained in the running Auto Sequence[®] is displayed, each equipped with the progress bar and percentage and the test status information. Individual results with status are displayed in real time. Click on any single test name to expand its window and view the results during test.

Currently running single test name is displayed in the field at the top of the Results window. It also contains the flow control commands **S**.

- **2** Summary: statistics of the running test is presented in real time. All individual results of all single tests are included in the count.
- OPProgress bar and percentage: the total progress of the Auto Sequence[®] is calculated and displayed in the top right corner of the Commands area. Completion of each single test is displayed in its name row together with it status.
- • Notifications of the Auto Sequence[®] test flow appear at the bottom of the Results sub window. Some require user action on test flow control.
- **6** Commands for user test flow control:

0

- **Provide a step**: continue with the next single test after the user action is completed.
- Next test: immediately stop the currently running single test and jump to the next step.

> 5 5		RUNNING MEASUREMENTS 70% 🖗				
RESULTS	0	APPLIED PART LEAKAGE	(DIRECT)			
 Visual IEC/EN 62353 	•		•			
 Protective Earth Resistance 			· 100%) •			
 Protective Earth Resistance 			· 100%) •			
 Protective Earth Resistance 						
 Applied Part Leakage (direct) 			50%			
~ I vMn LX1			0.001 mA 🌘			
✓ I VMn LX4			0.000 mA 🔹			
~ I VMn LX7			0.001 mA 🔹			
~ I VMr LX1			no value 🔿			
∽ vMr LX4			no value 🜼			
~ I vMr LX7			no value 📀			
			THE DEVICE UNDER TEST!			
SUMMARY						
Elapsed time	00:01:09	Passes	7			
Started at	14:03:05	Fails	0			
Estimated time to finish	00:00:31	Remaining	3			
		Total	10			

Stop test: immediately stop the Auto Sequence[®] test.

Figure 9.6: Auto Sequence[®] test progress data and notifications

Early interrupt

When the running single test is interrupted with the **Next test** command in an early stage, before any individual result are acquired, the user is prompted to manually select the further flow of the Auto Sequence[®] test, see **Figure 9.7**. Manual flow control menu appear only when all *Operation after end of test* flow commands of interrupted step are set to Manual.




Figure 9.7: Manual Auto Sequence[®] flow control prompt

Optional decisions available are:

- Stop test: Auto Sequence® test is stopped.
- Next step: proceed with next Auto Sequence® test step.
- **C** Repeat step: the interrupted test is repeated.
- End loop: exit the test loop where you have multiple steps in the measurement and proceed with the next measurement.

9.2.4 Auto Sequence results

The test screen of a finished Auto Sequence[®] (**Figure 9.8**) provides the following information and options:

- • The overall Auto Sequence[®] test state and status.
- **2** Auto Sequence[®] test steps' completion percentage and status.
- B Results: individual results of the Auto Sequence[®] test steps, each with status. Click on the step name to open its individual results list and scroll up/down to survey all.
 - Summary: Auto Sequence® test statistics.
- **D** Option commands.

AP-BF_Vis_PE_Leak	>			FINISHED
PROPERTIES				
© CONNECTIONS	RESULTS			
×	~ Visual IEC/EN 62353			
Visual IEC/EN 62353	~ Protective Earth Resistan	nce	2	100%
	~ Protective Earth Resistan	nce	0	
Single test	^ Protective Earth Resistar	nce		
OPERATION AFTER END OF TEST				
▼	~ R			0.00 Ohm 🌘
Protective Earth Resistance				
Single test	 Applied Part Leakage (d 	irect)		(
Operation after end of test				
v	~ L VMn LX1		6	0.001 mA 🔹
Applied Part Leakage (direct)	~ I VMn LX4		•	0.000 mA 🔹
Single test	∼ I vMn LX7			0.000 mA 🔹
OPERATION AFTER END OF TEST	~ Lan. 194			0.000 mA 🔹
	SUMMARY			
Result	Elapsed time	00:01:01	Passes	
Result Result SCREEN	Started at	09:23:25		
NESULT SCREEN	Finished at	09:24:27	Remaining	
			Total	

Figure 9.8: Auto Sequence® test results screen

Options after the Auto Sequence® test execution is finished:

- Save the Auto Sequence[®] test results to the location in the data structure and close the window.
- Cancel : Discard results, the Auto Sequence[®] test results are not saved. The user is prompted to confirm the action by selecting the option from the confirmation window, **Figure 9.9**.



Figure 9.9: Discard Auto Sequence[®] test results confirmation

Retesting using Auto Sequence® tests

When the finished Auto Sequence[®] test is opened and executed again, the new results are saved at the same location as the previously finished Auto Sequence[®] test. Confirmation is necessary. Previous Auto Sequence[®] test results remain unchanged.

9.3 Standardized inspections

Standardized inspections can be used as single test measurements or as measurements step within an Auto Sequence[®] test.

9.3.1 Standardized inspection test in the user workspace

Standardized inspection test window can be opened using different methods:



Inspection New Inspection test: Click on the Inspection icon in the Measurement group of Structure tab to open a list of available standardized inspection tests. Select one by clicking on its name. The inspection window opens in the user workspace. Inspection test can be saved and/or executed on PC.



Run Run an inspection test: Select an inspection test appended to the data structure in the tree view and click the Run command in the inspection test line, or in the Testing menu tab. An opened inspection test can be executed on PC or saved.



• View :-Inspection test: Select an inspection test appended to the data structure in the tree view and click on the View command within inspection test in the inspection test line, or in the Testing menu tab. An opened Inspection test can be surveyed only.

9.3.2 Executing the standardized inspection tests

Standardized inspection test window is presented on Figure 9.10:

- Inspection test name.
- 2 Header with the inspection test overall status and status selection list on the right.
- 3 Inspection test items with status selection lists on the right.
- Utem status selection list. Click on the triangle to open it.
- **5** Option commands




Figure 9.10: Standardized Inspection test window

Executing the standardized inspection:

• Perform each item in the inspection test and manually select the established status from the list.

- Overall status of the inspection test is set automatically, derived from the statuses of the inspection Items.
- Overall status of the Inspection test can also be set manually.
- **OK command:** confirm changes and save the inspection test results.
- **Cancel:** discard the changes and close the inspection test window. User is prompted to confirm cancelation.

Overall status automatic set rules:

- Pass: At least one Inspection item has status Pass and there are no Fail or Empty item statuses.
- Fail: At least one Inspection item has status Fail, regardless of the rest of items statuses.
- Empty: At least one Inspection item has status Empty, regardless of the rest of items statuses and there are no Fail item statuses.
- Checked: All inspection items have status Checked.

10 Measurement Organizer

Measurement Organizer is a tool for managing the user defined Auto Sequence[®] tests and inspections. Auto Sequence[®] tests can be organized into dedicated groups (*.atmpx files), saved to the PC file folder and uploaded to the instrument. A custom inspection is stored as a file with the *.indf ending. An existing Auto Sequence[®] group file can be reopened or downloaded from the instrument for further editing.

The Measurements Organizer can be accessed from the Tools group menu of the Home tab.

10.1 Measurements Organizer main menu

The Measurements Organizer opens in a separate window. The Main menu tab is active by default, **Figure 10.1**. Functionality and user options are:

• **OPEN FILES:** File storage location selection:

	() Recent files
0	() necent mes
0	Computer
0	Cloud
0	Instrument

0

• **2** Navigation tool: menu adopted to the selected file storage location.

- Recent files: List of recently edited files.
- Computer: browse the PC folder locations for files to open.
- Cloud: browse Cloud storage for files to open.
- **Instrument:** browse for sequences saved on the instrument. Set the communication port and select the Browse command to access the list of files saved in the instrument memory. Select a file and click Edit command to open it in the Auto Sequence[®] editor.
- **3** Create: Auto Sequence[®] editor and Custom inspection editor tools.

Measurements Organizer				-	×
e					
Main menu					_
Main menu Welcome to Measurements C To get you started quickly, we	brganizer e have provided some links to the most important features	8			
	Cloud 2		Create		
() Recent files	← Location: /	C ¢	+ Auto Sequence® group Create new Auto Sequence® group		
	Hospital Equipment periodic test				
Computer	Intenzive care periodic tests		Custom inspection Create new Custom inspection definition		
	KB_MZ_100				
Cloud	Machine test PAT test				
Instrument	Production test				
	test				
	Hospital equipment safety				
	Mains installations safety				
	Portable appliance safety				

Figure 10.1: Measurement organizer Main menu tab

10.2 Auto Sequence[®] Editor Workspace

Select Auto Sequence[®] group from the Main menu tab. A new Auto Sequence[®] appears on the screen in an empty Editor workspace in a separate tab. All necessary tools for creating Auto Sequence[®] tests-are provided in the three main Editor areas, see **Figure 10.2**:

- Auto Sequence[®] group is displayed as a folder tree structure. By default, it only contains a single new Auto Sequence[®]. See **Chapter 10.2.1** for details.
- **O** Auto Sequence[•] content view: the new Auto Sequence[®] template contains a header field PROPERTIES and footer field RESULT SCREEN by default. User can add single tests, inspections and flow commands to complete the Auto Sequence[®] test. See **Chapter 10.2.1** for details.
- Single tests tab and inspections tab. Each contains the available single tests and inspection tests organized in groups. Click on the triangle in front of the group () to open the list. Click on the sign in front of the selected test to add it to the currently active Auto Sequence[®] test.
- Summary: the selected Auto Sequence[®] test statistics, estimated run time and the total count of individual tests contained in all of the single tests.
 - **5** File: access the commands for:
 - o Standard file managing commands
 - **Upload**: open the communication settings window for uploading the file to the instrument.

Main menu New Auto Sequence® ×				
 New Auto Sequence[®] 1 2 3 1 	New Auto Sequence® PROPERTIES RESULT SCREEN 2	ROW COMMAND X ROW COMMAND X	* 200 MULUEMENT	Single tests Inspections All measurements MEIEC/EN 60601 MEIEC/EN 62353 Portable appliances COMPLETE ECG IBP RESP Temperature Mains Voltage Touch Voltage Voltage P-P
			En	N. C

Figure 10.2: New Auto Sequence[®] editor start screen tab

10.2.1 Auto Sequence[®] group tree structure

See the image **Figure 10.3** for an example of a composed Auto Sequence[®] group, with elements organized into a tree structure.

The tree structure elements are folders containing Auto Sequence® tests.

Measurement Organizer




Figure 10.3: Auto Sequence[®] group tree structure

Commands for composing and editing an Auto Sequence[®] group tree structure are located at the left side, or in the line with the selected element:

- Add new folder: select a location in the tree structure and click to add a new subfolder
- Add new Auto Sequence[•]: select a location in the tree structure and click to add a new Auto Sequence[•]
- Create shortcut: select an Auto Sequence[®] and click to create its shortcut.*
- Delete selected element(s): click to delete the selected folder with all content or the selected Auto Sequence[®].
- o Edit: edit the element name. Confirm with click on ♥, cancel editing with ▼.
- Handle: click and drag the element to the desired location in the Auto Sequence[®] group tree structure. To move to the subtree locations first click on the triangle (*) in front of the folders.

*Note:

Auto Sequence[®] Shortcut is useful when same Auto Sequence[®] with the same settings should be located at two or more different folders in the tree structure to keep the organisation consistent. Shortcut after creation is located beside the original Auto Sequence[®]. Select the Auto Sequence[®] shortcut handle and drag it to the desired location.

When a shortcut is selected, contents of the original Auto Sequence[®] is displayed. All additional editing and settings changes of the original Auto Sequence[®] are reflected in its shortcuts too.

10.2.2 Auto Sequence[®] content view

The selected Auto Sequence[®] is presented in detail in the Content view, see **Figure 10.4**. It is the main user working area. All necessary editing fields and commands are provided to adapt the Auto Sequence[®] to its purpose and the object under test. An Auto Sequence[®] is presented as a series of mandatory and user-inserted steps. Each step is presented in a separate field:

- **Header**: mandatory beginning step that enables setting the general PROPERTIES of the Auto Sequence[®] and configuring CONNECTIONS of the tested equipment.
- **O** Measurement steps: user-inserted single test and inspection test steps.
- B Result screen: mandatory ending step for setting the actions at the end of the Auto Sequence[®] and viewing the results.

	PROPERTIES			
0	Connections			
	Ψ			
Visual IEC/EN 62353		+	FLOW COMMAND	0
	Single test			
Operatio	ON AFTER END OF TEST			
Protective Earth Resistance	è	+	ELOW COMMAND	0
0000000	SINGLE TEST			
OPERATIO	W			
Applied Part Leakage (direct)		+	FLOW COMMAND	¢
	Single test			
Operatio	ON AFTER END OF TEST			
	à	+	FLOW COMMAND	

Figure 10.4: Auto Sequence® Content view

10.2.2.1 Auto Sequence[®] Header step options

Header step options are located in:

- Heading row
- Properties field
- Connections field

Heading row options

- AP-BF_Vis_PE_Leak Auto Sequence® name at the left side. It can be renamed in the tree structure.
- Add a user selectable flow command to be executed at beginning of the sequence.

Properties field setting options

Click the PROPERTIES field of the Header step to expand it and configure the options (Figure 10.5):

- Test optimization is checked as default to enable:
 - Individual medical leakage current tests will not be carried out one after another. They will be combined in a manner to achieve minimum power-ons and power-offs of the tested equipment and therefore the fastest run time.
 - The flow of the Auto Sequence[®] will change appropriately. Some options will become inactive to prevent any unnecessary delays and power-ons or power-off-s.

• A group of single test steps that are subjected to optimization is indicated with grey background and marked with the ¹⁰⁰ symbol, see **Figure 10.6**.

Note: It is recommended to have the test optimization always enabled.

		OPERTIES		
🕑 Test o	ptimization			
Auto Sequ U456	ence®code			
Descriptio Patient r				
	HR			
J-	~ 75		-0-	
1	BP 120/70		0	-0-

Figure 10.5: Properties field setting options with added image

- Auto Sequence[•] code: A user defined short code can be entered. The code is printed on the labels / tags and can be used as an identification of the Auto Sequence[•] in case of retesting.
- **Description:** A short description of the Auto Sequence[®] can be entered.
- Add Image: An image can be added.

Earth Leakage	🕇 FLOW COMMAND 🌣 🗙	: O Earth Leakage	+ FLOW COMMANE
Single t	EST	Single T	FST
OPERATION AFTER	END OF TEST	Operation after	
Pause		Pause	
Patient Leakage (to earth)	+ FLOW COMMAND & ×	Patient Leakage (to earth)	+ FLOW COMMAND
Pause		Pause	
Single t	EST	Single t	EST
Operation after	END OF TEST	OPERATION AFTER	END OF TEST
▼		· · · · · · · · · · · · · · · · · · ·	

Figure 10.6: Example of non-optimized (left) and optimized (right) test steps

Connections field setting options

The Connections field appears in the Header field when the first single test is added to the Auto Sequence[®] template. It provides the configuration for equipment under test, **Figure 10.7**. Click to the field name to access the subfields:

- PARAMETERS: Consider the test object's protection class (Class), presence of functional earthing (FE), signals inputs / outputs (SIO), the time it takes for the tested equipment to switch ON / OFF. Select the appropriate parameters from the lists. Note that the Auto Sequence[®] will not upload to the instrument or start unless all parameters are filled in.
- **CONNECTIONS:** Consider the test object's applied parts type, their functions and number of connections in each of them. Use commands to:
 - ADD CONNECTION: Select a connection type from the drop-down menu. Each connection type is presented in a sub field.
 - Expand the added connections to view additional information and functions.

- ADD FUNCTION (only applicable to type BF): Add a connection sub-function (SF). It is displayed in the type's subfield.
- ➤ Delete connection / function: Available on mouse-over the connection.
- Handle to reorder connections / functions; available on mouse-over the subfield name. Click and drag the selected connection to its desired location.
- **Rename** connection / function: select its default name and type in the custom name.

7	CONNECTIONS	
PARAMETERS		
Class		•
FE		Yes 🔻
SIO		Yes 🔻
Test On delay		5 s 🔻
Manual power off		Off ▼
Power Off delay		5 s ▼
CONNECTIONS		
CF AP CF	• +	
► CF AP CF	• • • +	
ADD CONNECTION		
AP AP CF		
AP AP CF		
	000000	
	4 5 6 7 8 9 10	

Figure 10.7: Connection setting options - default names (left) and custom names (right)

Note:

For more information about setting Connections refer to the instrument's user manual.

10.2.2.2 Auto Sequence® Measurement step

The Measurement step can be a single test measurement or an inspection test. Each is presented in a separate field where all the necessary settings and commands are provided to for adapting the measurement to it purpose in the Auto Sequence[®] test.

Measurement step options (Figure 10.8) are:

- Header row:
 - o AP-BF_Vis_PE_Leak Single test name at the left side.
 - Add a flow command field. See the Chapter 10.2.3 Flow commands for details.
 - Enter the Step count number to repeat the measurement more than once.
 - Remove the measurement step from the Auto Sequence[®] test.
- SINGLE TEST field: click to open single test subfields:
 - PARAMETERS: Set the single test parameters according to Auto Sequence[®] test requirements.
 - o LIMITS: Set the test limits according to Auto Sequence® test requirements.
 - o RESULTS: View the results and their details.

- **OPERATION AFTER END OF TEST field**: The user can select the action after the measurement is finished with respect to the result status.
 - o Auto: The Auto Sequence automatically moves to the next step.
 - Manual: Manually confirm moving to the next step, or use other commands during the pause.

	Single test	ADD FLOW COMMAND
✓ RESULTS		Buzzer mode
PARAMETERS		Inspection expert mode
Duration		No notifications mode
V mains		Pause
Condition		All 🔻
APs		All 🔻
NEPs		All 🔻
FE		All 🔻
Result		All 🔻
Test (AP=>PE)		All 🔻
Umax (calc)		Mains 🔻
Limit (CF,NC,ac)		10 uA 🔻
Limit (CF,NC,dc)		10 uA 🔻
Limit (CF,SFC,ac)		50 uA 🔻
Limit (CF,SFC,dc)		50 uA 🔻
Оре	RATION AFTER END OF	TEST
Operation after end of	test - Pass	Auto 🔻
Operation after end of	test - Fail	Auto 🔻


Figure 10.8: Measurement step setting options (Flow command selection field opened)

10.2.2.3 Auto Sequence[®] Result screen step

The Result screen is the mandatory final step of the Auto Sequence[®] test, Figure 10.9.

• • COMMAND : Add flow command field, see the Chapter 10.2.3 Flow commands for details.

	+ FLOW COMMAND X
RESULT SCREEN	ADD FLOW COMMAND
	Buzzer mode
	Inspection expert mode
	No notifications mode
	Pause



10.2.3 Flow commands

Flow commands are used to control the flow of measurements.

They are accessible via 🔸 ROW COMMAND , 🌣 , 🕨 buttons.

10.2.3.1 Number of measurement steps

Often the same measurement step must be performed on multiple points on the device under test. Set how many times a measurement step will be repeated. All steps are stored in the Auto Sequence® result as independent measurements.

Edit step properties	, 		×
Step count		1	

Figure 10.10: Measurement step repetition number setting

10.2.3.2 Pause

A Pause command with a text message or a picture can be inserted anywhere in the Autosequence[®]. Warning icon can be used alone or added to a text message. Custom text message can be entered into the field Screen Text, see **Figure 10.11: Pause setting examples**.

Parameters:

- Pause type:
 - Show text and/or warning. Enter text and add the warning icon (optional).
 - o Show picture: selected picture will be shown. Image path: Browse to location.
- Duration: Number in seconds or infinite.

			Pause
	Pause	Pause type	ShowTextAndOrWarning
Pause type	ShowPicture 🔻	Duration	5
Duration	5	Screen text	demo text
Image path		Show warning icon	0

Figure 10.11: Pause setting examples

10.2.3.3 Buzzer mode

In Buzzer mode State = On, passed or failed measurement is indicated with beeps.

- Pass double beep after the test
- Fail long beep after the test

Beep happens right after the single test measurement is finished.

10.2.3.4 No notifications mode

No notification mode setting influence pre-test warnings presentation.

- State = On, the pre-test warnings will be skipped.
- State = Off, the pre-test warnings are presented.

10.2.3.5 Inspection Expert mode

This flow command enables automatic application of inspection tickers.

- State = On the Visual inspection screen and Functional inspection screen within Auto Sequence[®] are displayed for 1 second and an overall PASS is automatically applied at the end of inspection test.
- State = Off automatic applying of tickers is disabled.

10.2.3.6 Operation after end of test

This flow command controls the process of the Auto Sequence[®] in regard to the measurement results. The operation can be individually set for the case the measurement passed, failed or ended without a status.

- Manual: The test sequence stops and waits for appropriate command (RUN key, external command...) to proceed.
- Auto: The test sequence automatically proceeds

10.3 Custom Inspection Editor Workspace

Create a custom inspection by choosing the Create Custom Inspection command in the Main menu of the Measurement Organiser. A new inspection with the default contents opens in the editor workspace, see **Figure 10.12**.

Main menu 🗦 New Cu	stom inspection ×		
	 Custom Inspection 		
	Custom Inspection	Pass_Fail_Checked_Empty ▼	

Figure 10.12: Custom inspection Editor Start screen

Use the *Save as* command from the Measurement organizer File menu (**Figure 10.13**) to rename and save the Custom inspection file to the selected location. File menu can also be used to open an existing file.

File		
	New	,
	Open	Ctrl+O
	Save	Ctrl+S
	Save As	
	Save to cloud	
	Recent files	•
	Close	Ctrl+W
	Exit	Ctrl+X

Figure 10.13: Measurement organizer File menu

Inspection field contains a Header at the top and a single inspection item as the Header child. All necessary tools are provided to configure the Custom inspection to it purpose.

Tools available within Inspection field Header are:

- Custom Inspection == Inspection name: Select it and type Custom inspection name.
- Add Child: click to add an inspection item.
- Remove: click to remove an inspection step;
 - Undo prompt appears on the screen for a short time. Click on it to cancel the inspection step removal.
- **Open / Collapse** the inspection content: click on the arrow on the left of the inspection name.
- **Handle**: click handle in the middle of the header and drag it to rearrange the inspection order. The handle appears only when the file contains multiple Inspections.

Tools available in the Inspection Item row are:

- Custom Inspection Item name: Select it and type in the custom item name.
- Add Child: click to add a child item to the selected item. It creates a tree-like structure.
 Pass_Fail_Checked_Empty
 × ||

```
Pass_Fail_Checked_Empty
```

- **Type** of the item status: click on the arrow on the right side to open the item status selection options: *Pass_Fail_Checked_Empty* or *Pass_Fail_Empty*.
- **Remove:** click to remove an inspection item.
 - o Undo prompt appears on the screen for a short time. Click on it to cancel the inspection item removal.
- Handle: click and drag an item to the desired location in the Inspection field.

Tool available under the inspection field:

• **Add Inspection:** click to add an inspection. Custom created inspection file can contain multiple inspections. See **Figure 10.14**. for an example.

Measurements Organ	lizer			(S 1 _2)	
le					
🕈 Main menu 🛛 🎘 🕻	Custom Combo Inspection VS monitor ×				
	 Vital sign monitor visual 		×		
	 Mains supply 	Pass_Fail_Checked_Empty 🔻			
	no defective or damaged cable	Pass_Fail_Empty 🔻			
	no defective plug or plug pins	Pass_Fail_Empty 🔻			
	no defective or damaged mains switch	Pass_Fail_Empty 🔻			
	▼ Enclosure	Pass_Fail_Checked_Empty 🔻			
	no missing or damaged parts	Pass_Fail_Empty 🔻			
	no defective accessory connectors	Pass_Fail_Empty 🔻			
	correct marking and warnings	Pass_Fail_Empty 🔻			
	▼ LCD screen	Pass_Fail_Checked_Empty 🔻			
	no cracks or scratches	Pass_Fail_Empty 🔻			
	 Vital sign monitor functional 		×		
	 Mains supply 	Pass_Fail_Checked_Empty 🔻			
	Correct Mains switch operation	Pass_Fail_Empty 💌			
	Correct Mains ON signalling	Pass_Fail_Empty 🔻			
	▼ LCD screen	Pass_Fail_Checked_Empty ▼			
	no vertical lines or dead or bright pixles	Pass_Fail_Empty ▼			
	no flickering or jittering	Pass_Fail_Empty 🔻			

Figure 10.14: Custom Inspection file with multiple inspections Editor Workspace

Note:

Custom created inspection tests cannot be used as single test measurement appended to the test structure element. They can only be used as a step in the Autosequence[®] test.

11 Printing results

Metrel Medical ES Manager provides a simple print format for the test project results printout. The command Print results is located in the Tools group of the Home tab.

11.1 Printing results options

Create a header of test project results printout form:

- Open the Settings menu located in the Main tab menu.
- Select the Report settings tab.
- Printout results tab (Figure 11.1) header setup options are:
 - **Change logo image**: Opens the file browser to navigate to the logos' image folder. Select the image and click the Open command.
 - Click **OK** command to confirm Logo change or **Cancel** to discard changes.
 - Clear logo image: Delete the existing logo. Command cannot be undone.

Settings									
General	Language	Report settings	Cloud						
Printout re									
Change	logo image		o image	l	_	0	g	0	
OK	Cancel								

Figure 11.1: Print results Logo setup

11.2 Printing results format

Procedure to open Print format form:

- Open the test project data structure file. In the tree view select the parent object that contains all the elements whose testing results are to be printed. Select the root node to print all results. Select a single measurement to print only its results.
- Click on Print Results in the Tools group of the Home tab to open the test project results document in the Preview manager window,
- Figure 11.2.

Print	K 🖽 - 🅀 Find	Previous Page	▶ □□ Q		-		
Quick Print	Scale B - B Thumbnails	First Next Page	📲 🗆 q.		- P	-	
] Parameters	- Editing Field		Many Pages	Page Color	· -	Close	
Print	Page Se 15	Navigation	Zoom	Page Background	Ex	Close	
		5		5 5			
	_						
	N	lew I		-			
			00				
			_0q				
		•					
		•					
	Instrument data:	0					
	Model: MediTest MI 6601	Serial Number: 2	0211663 Calibration	n date: 20/05/2020			
	User:		No. (all salates salates		_		
	PM01001 Demo PM Comme	ent: Class II Patient monitor	SenseVis IC Room 20/				
	AP-BF_Vis_Iso_Leak			Pass			
	Visual IEC/EN 62353			Pass			
	Riso LN-NEP			Pass			
	Riso LN-NEP	4		Pass			
	Riso LN-NEP			Pass			
	Riso LN-NEP			Pass			
	Riso LN-NEP			Pass			
	Applied Part Leakage (d	irect)		Pass			
	PM01002 Demo PM		SenseVis IC Room 20/	04/20 20/10/20 Pass			
	AP-CF_Vis_PE_Leak			Pass			
	Protective Earth Resista			Dees			
				Pass			
	Applied Part Leakage (d PM02001 Demo PM	irect)	SenseVis IC Room 20/				
			Sensevis IC Room 20/		/		
	AP-BF_Vis_Iso_Leak	^		Empty			
	Visual IEC/EN 62353	4		Empty			
	Riso LN-NEP			Empty			
	Applied Part Leakage (d	irect)	0.11.10.0	Empty			
	PM02002 Demo PM		SenseVis IC Room 20/		Y		
	AP-CF_Vis_PE_Leak			Empty			
				Empty			
	Visual IEC/EN 62353						
	Protective Earth Resista			Empty			
				Empty Empty			
	Protective Earth Resista						
	Protective Earth Resista						
	Protective Earth Resista Applied Part Leakage (d	irect)	5/0023				
	Protective Earth Resista		5/2023		1/1		
	Protective Earth Resista Applied Part Leakage (d	irect)	5/2023		1/1		

Figure 11.2: Print Results document format preview

- **Header** appears on the top of the first page:
 - Logo (**0**) is set in the Report settings tab of the Settings menu.
- Instrument data (2) gets set automatically when the data structure file is downloaded.
- Footer of the results printout (3) appears on the bottom of each page and contains:
 - o Signature field enter by hand
 - o Created date is set automatically
 - o Page numbering is set automatically.
- Measurement results are presented in table (④), each table starts with Appliance name, it path within data structure tree and overall status indication, followed with individual test result status.
- Edit the results page layout and print it using the preview manager commands.

Preview manager also provides the Export group with options:

- Export to PDF file: Create file to save document in PDF format.
- Export to Image file: Create file to save document in image file format (selectable from menu).
- E-Mail as PDF file: Create file in PDF format and attach it to E-Mail message.
- E-Mail as Image file: Create file in image file format and attach it to E-Mail message.

12Creating Reports

User can create, manage and archive Reports of Test project measurements with WebReports application. The application is a part of the Metrel Ecosystem platform, using the same user repository as Metrel Cloud. Before first use, new user must first register Metrel Cloud account and then enter WebReports License key to Cloud account.

For an easy start different pre-defined report templates are offered by Metrel. They can be used as they are or be modified and used as personal templates. The tool enables to create customized reports suited to one's personal demands. New templates can be made from scratch.

WebReports templates are accessed from Metrel Medical ES Manager Home tab Reports menu:



- Web Reports : Web Report command is active when Data structure file with appended measurements is opened.
 - **Option**: Data structure file was opened from storage location other than Cloud: browser is opened first to navigate and save Data structure file to Cloud folder.
 - Option: Data structure file was opened from Cloud storage: Report Selector window with available list of graphic Report template files is opened, see Figure 12.1. A convenient tool is available to search for pre-made Report template.

Name	Date modified V	ersion
•Ds		
= Equals	21/04/2022	
(Deve est servel	25/03/2022	
	28/06/2022	
Contains	24/06/2022	-
• Does not contain	24/11/2021	
	22/11/2021	
stác Is like	22/11/2021	
** Is not like	23/03/2022 24/11/2021	
Rec Regine with	24/11/2021	
ouc Begins with	15/05/2023	
Ends with	15/05/2023	
> Is greater than	15/05/2023	
5	15/05/2023	
	26/11/2021	
< Is less than	30/03/2022	1
	15/05/2023	· · · · ·
S is less than or equal to	15/05/2023	
Test patient monitor 03	15/05/2023	
Test patient monitor 04	15/05/2023	
Test_patient monitor_05	15/05/2023	
Test_patient monitor_06	15/05/2023	
Test_patient monitor_07	15/05/2023	
z_strukturo	28/06/2022	

Figure 12.1: Report Selector window with Find tool opened

• Create report : Select appropriate Report graphic template and confirm. WebReports editor is opened.

Refer to WebReports User Manual for details on how to create templates and edit, manage and archive reports.

13Troubleshooting

13.1 Introduction to troubleshooting

This chapter contains information that might be helpful if you encounter problems while working with Metrel Medical ES Manager. If you don't find the help required to troubleshoot your problems, you can contact us directly. Check the <u>Online support information</u> for information on how to reach us.

Your comments will help us improve our product and update the manuals with any necessary information.

13.2 Instrument connection troubleshooting



Figure 13.1: Instrument connection error message

If the instrument settings are configured correctly, and the communication still does not work, please check the following:

- Make sure that the cable is not broken and that the connectors are inserted properly. If you have another device which can be connected using the same cable, consider testing the cable using that device to make sure that it works.
- Make sure that you close all other applications that might be using the port. You can also try restarting Windows to make sure that all such applications are closed.
- If none of this helps, there might be a problem with the device driver for the port you are using. You could try to reinstall or repair Metrel Medical ES Manager to make sure that your USB drivers are up to date.

13.3 Reporting problems to Metrel

This section deals with unexpected problems, errors and exceptions that may occur during usage of Metrel Medical ES Manager. If you have encountered an error while using Metrel Medical ES Manager, the best thing to do is to report this issue to our support at <a href="mailto: This will help us find a solution quickly, and provide you with an updated version, hotfix, or help you extract the data from your instrument if the issue prevents you from doing so.

13.3.1 Bug report checklist

This section contains a brief list of items which your bug report should contain. For detailed information about each step, consult the next section.

- 1. Regardless of the problem you are reporting, the minimum amount of information we need is:
 - A short description of the problem and, most importantly, steps needed to reproduce it;
 - Metrel Medical ES Manager log file (located in ...\Metrel\Medical MESM\logs.

140nline support information

14.1 How to get support?

14.1.1 Contact your local distributor!

Most issues can be solved by your local distributor. This support is accessible in your local language. Find the closest distributor here: <u>https://www.metrel.si/en/locations/</u>

14.1.2 Use the online contact form

You can also contact us using our <u>Contact form</u> [https://www.metrel.si/en/contact/]. That web page also contains a map describing the route to our Company's headquarters.

Appendix A Using Auto Sequence®s to test according to IEC/EN 60601

Metrel prepared the 4 Auto Sequences for testing according to IEC/EN 60601 with the intention of enabling the users to edit them to the requirements of any tested device as easily as possible.

In their default state, they contain every possible measurement mentioned in IEC/EN 60601 applicable to the **tested device's protection class. The actual choice of tests depends on the manufacturer's requirements and is** usually much shorter.

o ×

The following chapter deals with features of the default sequences, their editing, and the recommended workflow.

Getting started

Opening the default IEC/EN 60601 Auto Sequences group .atmpx file shows the window below.

menu HEC 60601 full length ×				
🗰 Class I, fixed supply wire	Cass I, fixed supply wire	· Martin and Assault	Single tests	Inspections
Class L removable supply wire	Properties		ADD MERCUREMENTS + All measure	rments
Class II Class IP	O CONNECTIONS		 ME IEC/EN 	
	*		► ME IEC/EN	
	Protective Earth Resistance		 Portable ap 	
	Siville YEST		COMPLE	TE
	OPERATION AFTER END OF	TEST	IBP	
			RESP	
	IO Earth Leakage	+ n.m. (1	Temperature	ure
	Sivale TET		Mains Vo	
	Oversation artists pilo or	1000	Touch Vo	ltage.
		0 x	Voltage P	je P-P
	Touch Current	0 ×		
	Single yest			
	Touch Current (NEP to NEP)	e x		
	Sindué rest			
	Operation after the of			
	Patient Leakage (to earth)	0 x		
	Single 1117			
	OPERATION AFTER END OF	1942		
	Patient Leakage (Vext on AP)	0 x		
	Sivale 7857			
	Operation after end of	1817		
	Patient Leakage (Vext on SIO)	0 x		
	Single test			
	OPERATION AFTER END OF	TELT		
		0 X		
	Patient Leakage (Vext on NEP)			
	Single test		SUMMARY	
	Overation artist too or			
	Total Patient Leakage (to earth)	0 <u>x</u>	Estimated time	Number of tests

Figure 2: Open Auto Sequence Class I, fixed supply wire

The Auto Sequences in the group are collected in the list on the left. The contents of the currently open Auto Sequence are in the middle. The right column contains a list of supported measurements that can be added to the Auto Sequences.

Metrel prepared 4 Auto Sequences. The choice depends on the protection class and supply wire:

- Class I, fixed wire: contains the Protective earth measurement with test current of 25A and limit 200 $m\Omega$. The measurement is done between any earthed conductive part and the ground contact at end of the supply cable. Uses the worst-case mains voltage to measure leakage currents. Add steps to measure additional earthed accessible parts using a probe.
- Class I, removable wire: contains Protective earth measurement with test current of 25A and limit of 100 mΩ. The measurement is done between any earthed conductive part and the inlet for cable. Uses the worst-case mains voltage to measure leakage currents. Add steps to this measurement in case of any earthed accessible parts.

- Class II: **Doesn't contain the Protective earth measurement, uses the worst**-case mains voltage to measure leakage currents.
- Class IP: Doesn't contain the Protective earth measurement. Uses the worst-case supply voltage to measure the leakages. Set the Umax(calc) in all leakage current parameters to the supply voltage.

The first step of the Auto Sequence lists the Properties of the tested device. They must be filled in before testing according to the tested device and connections to the MI 6601.

- 1. Auto Sequence name
- 2. Add a Flow command. See the Chapter 10.2.3 Flow commands for more information on Flow commands.
- 3. **Test optimisation: allow MESMM to change the result sequence to minimise tested unit's resets and** streamline the testing. See MESMM user manual for more information.
- 4. Insert the desired Auto Sequence code. The code is printed on labels / tags and can be used as identification of the Auto Sequence[®] in case of retesting.
- 5. Tested device description.
- 6. Add image, often used for help screens or schematics.
- 7. Connections between the tested unit and the instrument. MESMM will automatically setup the necessary steps and limits for the leakage current tests. Any tests that are not applicable to the selected connection configuration will be automatically skipped.
- 8. Set protection class of the device.
 - a. Class I protected by single layer insulation and protective earthing.
 - b. Class II Protected by double insulation.
 - c. Class IP IP stands for internally powered. Covers all devices with batteries and special chargers.
- 9. Set presence of FE (functional earth) and SIO (Signal Input/Output) connections.
- 10. Set Test on delay. We recommend measuring the time the instrument needs from pressing the On button to full functionality, plus a couple of seconds to reach the On button.
- 11. Set Manual power off and Power off delay. If the Manual power off is off, this is the exact time the tested device takes to shut down. If the Manual power off delay is on, add some seconds for pressing the Off button manually.
- 12. Connections setup. Click the plus symbol to add the connections to the MI 6601 applied parts connectors. The possible uses of these connections are applied parts (AP), earthed parts (EP) and non-earthed parts (NEP). Set them in the same way as you connect the tested device to the MI 6601.

lass I, fixed supply wire	FLOW COMMAND
	Properties
Test optimization	
Auto Sequence®code	
Description Device with protection of	class I and fixed supply wire.
🕂 ADD IMAGE	
Parameters(s) with undefin	© Connections
PARAMETERS	
Class	1 🔻
FE	Yes 🔻
SIO	Yes 🔻
Test On delay	not set 🔻
Manual power off	Off 🔻
Power Off delay	5 s 🔻
CONNECTIONS	
+ ADD CONNECTION	
0.00	0000000
	3 4 5 6 7 8 9 10

Figure 3: Connections

Measurements are automatically populated with the required results and limits. If there are no results to measure, a warning will come on. User action is necessary to skip the measurement.

Measurements and parameters

See chapter 10.2.2 Auto Sequence[®] content view for more details on Auto Sequence content manipulation.

The sequence is designed so that the user can easily remove the unnecessary measurements and would only need to add any in case measurements by other standards are needed.

The manufacturers typically require only some of the leakage measurements available in the sequence. The rest can be simply removed using the X symbol in the upper right-hand corner. Save the Auto Sequence under a new name. The Metrel Auto Sequences are locked and cannot be overwritten.

Next to it is the Settings button, where the number of repetitions of the same measurement can be set. For example, set it for Protective earth resistance in case there are multiple earthed parts on the housing of the tested device. Consider that it will be necessary to use a probe to connect the each of the measured points.

The default state of the parameters in the leakage current measurements is to include the largest possible number of measuring points. Duration is set to the shortest possible time, 2s. Set the parameters as required by the tested device manufacturer. If uncertain, leave as is.

Example: patient monitor Philips MP70

The monitor has 10 connections for ECG that need to be tested, a functional earth, and an SIO connector that is grounded and therefore irrelevant for leakage current testing. The supply cable is fixed. The connectors are set into an earthed conductive board that can be covered by a lid, but the clasp is easy to detach, so we count it as accessible. There are 2 accessible screws on the back of the housing. There are no non-earthed conductive parts.

Setup of the Auto Sequence

The Auto Sequence chosen is the Class I, fixed wire. The connections window includes presence of FE, but not SIO. The applied parts can be considered a single part with 10 connections. Test on delay and power off delay were measured to 20 s and 5 s respectively.

ilips InteliVue MP70		+ FLOW COMMAND
	Properties	
	© Connections	
PARAMETERS		
Class		•
FE		Yes 🔻
SIO		No 🔻
Test On delay		20 s 🔻
Manual power off		Off ▼
Power Off delay		5 s 🔻
CONNECTIONS		
► CF AP		
+ ADD CONNECTION		
AP		
		10

Figure 4: Example of connections for MP70

This setup is enough for a thorough testing. However, requirements are less extensive and time can be saved by shortening the measurement list.

Protective earth resistance stays, it is the most crucial test in Class I device. There are 3 measuring points, so the steps are set to 3.

AP	F. P. Handard and State		
•	Edit step properties	×	
1	Step count	3	
Protective Earth Res		ND 🗘 🗙	
O Patient Leakage (Vext	on AP)	🕈 FLOW COMMAND 🌣 🗙	

Figure 5: Inputting the step number

Earth leakage, patient leakage to earth and total patient leakage to earth are not required by the manufacturer in this case, as referenced in the Service manual. The SIO is irrelevant, and there are no NEPs, which leaves a lot of the remaining measurements empty. This is the cleaned-out sequence on Figure 5: Inputting the step number.

Philips InteliVue MP70	+ FLOW COMMAND ×
PROPERTIES	
CONNECTIONS	
Protective Earth Resistance	🛉 FLOW COMMAND 🌣 🗙
Single test	
Operation after end of test	
Ψ	
© Patient Leakage (Vext on AP)	🕈 FLOW COMMAND 🌣 🗙
Single test	
Operation after end of test	
Total Patient Leakage (Vext on AP)	🕂 FLOW COMMAND 🌣 🗙
Single test	
OPERATION AFTER END OF TEST	
Ψ	
	FLOW COMMAND X
RESULT SCREEN	